

PHILIPS

AVENT

Handle for breast
pump

Natural

Comfort breast pump

Manual

CP9282



To make the breast pump diaphragm move.

Find matching products on the specifications tab

Pushing this handle up and down creates suction. You can control the suction by varying the pressure on the handle.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCF330/40, SCF330/70, SCF330/20, SCF330/71, SCF330

