

PHILIPS

Adjustable crank set

ReActiv



PTE4010CA



Upgrades bikes to fit all users

Best for folks with limited mobility

The Adjustable crank set makes exercises in smaller range of motion possible on bikes. This accessory gives greater safety and accessibility to users with impaired range of motion as well as helping to reduce fear of pain.

Highly adjustable for different range of motion

- Crank length is adjustable for limited mobility

Highlights

Adjustable crank set

For users with limited mobility, the optional adjustable cranks are indexed for repeatable placement and have a starting range of motion

as small as 15 degrees. As their performance improves, they can continue their recovery to full range of motion.

Specifications

Accessory set

set: 1

Thread size

inches: Standard 9/16" thread

Adjustment degrees

From 15 degrees to full range: -

Applicable products

Upright bike

Recumbent bike

Pedal position adjustments

10 incremental: -

