

# PHILIPS

Full Face Mask (M)

DreamWear

Full Face Pack  
with Medium cushion  
and Small and Medium Frame



HH1126



## Like wearing nothing at all.\*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

#### **Magnetic clips**

- The headgear clips and mask cushion contain magnets

#### **Open field of vision**

- Open design offers a wide field of vision

#### **Innovative design**

- Freedom of movement
- A modular design for more options

#### **Comfort and ease of use**

- Easy to use
- Soft, flexible frame
- Under the nose full face cushion

# Highlights

## Freedom of movement



By moving the tube location to the top of the head, DreamWear Full Face offers more freedom of movement during the night. Users stated they didn't need to choose a sleeping position when wearing DreamWear Full Face.\*

## No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.\*\*

## Unobstructed view



Open design makes it easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.\*\*

## Like wearing nothing at all



The soft, flexible frame is comfortable to wear at night. Users have stated it feels like wearing nothing at all.\*

## Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

## Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

## Magnetic clips



Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

# Specifications

## Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

## Intent of use

To be used on patients <gt/>30 kg: with CPAP prescription

## Cleaning: Headgear, tubing

Hand wash with mild detergent: Rinse thoroughly, hang to dry<sup>1</sup>

## Cleaning: Mask

Hand wash with mild detergent: Rinse thoroughly, air dry<sup>1</sup>

## Replacement

Inspect daily for wear: Replace at first sign of wear

## Product details

Materials: Nylon, Polyurethane foam, Polycarbonate, Polyester, Silicone, Spandex  
Warranty: 90 days

Operating pressure: 4 to 30 cmH<sub>2</sub>O

Sound level: 17 dBA

Number of parts: 5

