



PHILIPS

Full Face Magnet
Clips (2 pack)

DreamWear

Full Face Replacement Part
1 set of magnetic clips (2 pack)



Product Picture

HH1132

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Easy to use

- Magnet clips for headgear attachment

Highlights

Easy to use



The magnet clips make attaching the headgear to the frame easy and fast!

Magnetic clips



Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Specifications

Warning: Magnetic interference

Beware with pacemakers: Read IFU's before use.

Product details

Warranty: 90 days

