

PHILIPS

Replacement
Mouthpiece Assembly

InnoSpire Go

Vibrating mesh technology

Replacement part

For use with InnoSpire Go



HH1345



Clinically proven vibrating mesh technology

The replacement mouthpiece assembly is designed to maintain the performance of InnoSpire Go and to efficiently deliver medication for children and adults.

Clinically proven vibrating mesh technology

- Vibrating mesh technology, used in hospitals worldwide

Easy to clean and maintain

- Only one part to clean

Provides 1 years use

- Provides 1 years use

Highlights

Easy to clean and maintain



The easily detachable mouthpiece makes the InnoSpire Go simple to use, clean and maintain.

Clinically proven technology



It consistently delivers aerosol in the optimal particle size, to target the medication to your lungs.

Therapy use



Giving you consistent and reliable therapy.

Specifications

Service

Mouthpiece assembly warranty: 1 year, recommend yearly replacement

Product details

Mouthpiece assembly: for use with InnoSpire Go

Maximum fill level: 8 ml

Mesh material: Nickel palladium

Outer body material: Grilamid™ (plastic)

Packaging size: 65mm x 65mm x 70mm

Packaged weight: 45g

