

PHILIPS

Fabric Wraps (2 pack)

DreamWear

Replacement Part

1 set (2 pack)

2 per pack



HH1120



Like wearing nothing at all.*

With the freedom and comfort of DreamWear, you can sleep how you want.

Easy to use

- Comfort and soft

Highlights

Soft fabric wraps



Provides a soft touch on your cheek

Specifications

Cleaning: Fabric wraps

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Polyester, Spandex

Warranty: 90 days



* * 2017 Philips User Preference Questionnaire

* * Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respiration Amara View and F&P Simplus; based on 10 days of use.

* ¹Please refer to Instructions For Use for more detailed cleaning instructions.