

Philips Sonicare Power  
Flosser 3000  
Cordless water flosser -  
Black

Removes up to 99.9% of plaque

2 flossing modes, 3 intensities  
Complete clean in 60 seconds  
Fast universal charging with USB-  
A



HX3826/33

## Flossing reinvented with QuadStream technology

Removes up to 99.9% of plaque\* in treated areas

Effortlessly thorough cleaning between teeth to improve your gum health. Unique Quad Stream tip cleans more surface area with less effort for more effective flossing every time, so it's easy to get flossing right.

### Easy and convenient

- Effortless technique
- Fast universal charging with USB-A to Small Plug cable
- Easy-to-fill 250 ml reservoir
- Standard tip for spot cleaning

### Faster, more effective flossing

- Complete clean in 60 seconds
- 2 flossing modes, 3 intensities
- Pulse Wave technology guides from tooth to tooth
- Quad Stream technology for faster, more effective flossing

**PHILIPS**  
**sonicare**

## Cordless water flosser - Black

Removes up to 99.9% of plaque 2 flossing modes, 3 intensities, Complete clean in 60 seconds, Fast universal charging with USB-A


# Highlights

### Quad Stream technology



Unique X-shaped Quad Stream tip separates flow into 4 water jets, which cover more surface area between teeth and along the gum line. Get a faster, deeper clean without the hassle of string floss.

### Pulse Wave technology



Gentle pulses of water keep you going and guide you from tooth to tooth in Deep Clean mode, so you get it right every time.

### 2 modes, 3 intensities



Choose the clean that suits your smile. Clean mode delivers continuous flow for effective everyday cleaning. Deep Clean pulses for more

thorough cleaning. Adjustable intensity level adapts for comfort.

### Complete clean in 60 seconds



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation. Get a thorough 360 clean in just 60 seconds from start to finish.

### Fast charging with USB-A to Small Plug cable



Convenient charging with USB-A to Small Plug cable. You can clean for up to 14 days on just one charge.


HX3826/33

### 250 ml reservoir



Reservoir capacity of 250 ml holds enough water for a recommended 60 second clean. Simply twist and remove or use the side-fill door for easy refilling.

### Effortless technique



Just point, press and clean for more effective flossing. Aim the tip between your teeth and along the gumline, and let it do all the work for you!

### Standard tip for spot cleaning



A standard single-stream tip boosts pressure for spot cleaning and removing food debris when you need it.

## Cordless water flosser - Black

Removes up to 99.9% of plaque 2 flossing modes, 3 intensities, Complete clean in 60 seconds, Fast universal charging with USB-A

HX3826/33

# Specifications

### Power

- Voltage: Multi-voltage charger

### Service

- Warranty: 2-year limited warranty

### Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

### Performance

- Cleaning: Clean entire mouth in 60–90 sec

### Items included

- Cordless Power Flosser: 1
- F1 Standard nozzle: 1

- F3 Quad Stream nozzle: 1
- Travel bag: 1
- USB charging cable: 1
- USB wall adapter: 1

### Modes

- Clean: For exceptional everyday cleaning
- Deep Clean+: For an invigorating deep clean
- Intensities: 3

### Health benefits

- Plaque removal: Removes up to 99% of plaque in treated areas\*

### Design and finishing

- Colour: Black



Issue date 2024-03-24

Version: 3.2.3

EAN: 87 20689 01113 6

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* In an in-vitro study. Actual results may vary.