

PHILIPS

Headgear
(Standard)

DreamWear

Standard



HH1119



Like wearing nothing at all*

With the freedom and comfort of DreamWear, you can sleep how you want.

Easy to use

- Soft and comfortable

Highlights

Innovative design

One modular soft fabric headgear for all frame sizes

Specifications

Cleaning: Headgear, tubing

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Nylon, Polyurethane foam, Spandex

Warranty: 90 days

