

Philips  
Massage Back Brace

VibWave deep shoulder massage

Constant hot compress  
3D elastic shoulder muscle aid  
Poor posture reminder

PPM4361



## You personal masseur and posture coach

### Healthy life with deep massage and health posture

This massager is light weight and wearable. You can enjoy deep vibration massage and hot compress to your shoulder at anytime. The smart built-in posture sensor can provide a kindly reminder when bad posture is detected.

#### Healthy materials applied

- Anti-bacterial and splash-proof fabrics

#### Deep relaxation for should muscles

- Constant temperature hot compress
- 3-area VibWave vibration massage for deep mucle relaxation

#### Assist for a better posture with back straight

- 3DSuit stereo-cut back brace assists should muscles
- Smart posture sensor to more easily improve posture habits

# PHILIPS

## Massage Back Brace

VibWave deep shoulder massage Constant hot compress, 3D elastic shoulder muscle aid, Poor posture reminder

# Specifications

### General specifications

- Product Dimension (LxWxH): 454x367x20 mm
- Net weight: 330 g
- Color: Blue

### Main parameters

- Working voltage: 3.7 V
- Working power: 12 W
- Battery: 5000 mAh

- Charging current: 2 A
- Charging voltage: 5 V

### Country of origin

- Made in: China

### Service

- Warranty period: 2 years

PPM4361DB/97

# Highlights

### Thermostatically hot compress

With temperature sensors built-in. The well controlled massage module provides a relaxing hot compress to improve the muscle blood circulation.

### 3DSuit assists should muscles

The stereo cut body fits your back well and the elastic strips help to reduce the burden of your shoulder muscles and straighten your back.

### 3-area VibWave vibration massage

The VibWave massage uses a high frequency vibration which goes deeply into muscle and relax muscle tissues efficiently.

### Smart posture sensor

The built-in posture sensor can work a whole day, providing a kindly reminder when bowed back is detected. Users will gradually get used to a healthier habit to keep back straight when sitting and standing.

### Anti-bacterial and splash-proof fabrics

The device is made of fabrics processed with nanotechnology. The anti-bacterial ingredients prevents the device to get smelly and unhealthy after use. The splash-proof layer keeps most dirt away and makes it easy to clean.



Issue date 2023-09-28

Version: 1.1.1

EAN: 48 97126 05053 0

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)