

Philips  
Knee massager

### Combined massage

Airbag full wrap  
wearable use in wireless

PPM5521



## Enjoy the moment to relax

### AirWave pressing and hot compress

Help relieve knee discomfort through air pressure and heat compress functions. A variety of massage modes and gears are available, suitable for different groups of people for knee massage and relaxation.

#### Combined massage, relieve discomfort

- Massage both knees together to relax
- 40-50°C hot compress temperature, 3 adjustable temperatures
- 4 massage modes, smart and convenient

#### Airbag full wrap design

- 10 airbag areas, air wave airbag massage
- Suspended three-dimensional support for the knee joint

#### Easy wireless use

- 30-minute intelligent timed massage, avoid excessive massage
- Built-in battery, charging use is available

# PHILIPS

## Knee massager

Combined massage Airbag full wrap, wearable use in wireless

# Specifications

### General specifications

- Color: Dark blue
- Net weight: about 0.78 kg
- Quantity of airbags: 10

### Main parameters

- Temperature of heating: 40-50 °C

- Battery: 5000 mAh
- Charging port: USB Type-C

### Country of origin

- Made in: China

### Service

- Warranty period: 2 years

PPM5521DB/97

# Highlights

### All-in-one knee massage

Dual massager design, both knees can enjoy knee care at the same time.

### Hot compress function

3 levels of hot compress temperature can be adjusted to promote blood circulation, effectively relieve muscle tension, reduce local pain, let people fully relax and relieve the fatigue of the day. Built-in precise NTC temperature control system, overheating automatic power off, to achieve intelligent anti-scald.

### Multiple massage modes

There are 4 preset massage modes, one button to start the massage, giving you the comfort you want.

### Airbag design

There are 10 three-dimensional airbags inside, which wrap the knees in multiple angles for massage, give the knee airbags a squeeze massage, and relax the knees

### FloatFit 3d knee pads

The knee part of the device is equipped with a ring-shaped air suspension, which fixes the patella by inflating the airbag and gives the patella a supporting effect. Patella massage can be achieved by switching between inflation and deflation cycles, relaxing the patella and relieving the pressure on the patella.

### Timing shut-off

A timed 30-minute working time for automatically shut off, effectively avoiding excessive massage and causing discomfort of redness and swelling of the skin.

### Use wireless

Built-in 5000mAh battery can be used for 120 minutes when fully charged. The device supports use in charging, so there is no worry about battery life.



Issue date 2023-09-28

Version: 1.1.1

EAN: 48 97126 05062 2

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)