


PHILIPS

AVENT

Single Electric Breast Pump
with Aftercare Support

With professional support

Includes 4 oz bottle

 Product Picture

SCF332



Philips are here to support your breastfeeding journey

With professional support

When you buy a Philips Avent single electric breast pump, you receive unlimited one-on-one video calls over 3 months with a WellVine breastfeeding professional to support you on your breastfeeding journey.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Easy to combine breast and bottle feeding

- Includes our Natural bottle and teat for natural latch-on

Other benefits

- Intuitive assembly. Easy visual matching of parts
- Compact lightweight design
- Compatible with other Philips Avent feeding products
- Easy cleaning due to the small number of separate parts
- Easy operation at the touch of a button
- Soft massage cushion with warm feel

Professional breastfeeding advice

- 3 months' unlimited video calls with breastfeeding experts

Highlights

Professional support

When you buy a Philips Avent electric breast pump, you get unlimited video calls over three months with WellVine Lactation Consultants to support you on your breastfeeding journey — wherever you are. WellVine Lactation Consultants are International Board Certified (IBCLC), the highest qualification in breastfeeding support. They have helped hundreds of mothers with advice on preparing for breastfeeding and what to expect when breastfeeding a newborn, as well as advice on breastfeeding issues, including low milk supply, breast pain, latching issues and more!

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means that you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Includes our Natural teat



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Soft massage cushion

Our massage cushion has a new soft, velvety texture that gives the skin a warm feeling for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let-down.

Compact design

The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. Easy to use on the go with batteries.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

Easy operation

A single electric breast pump is ideal for mums who pump milk frequently and value being able to pump effortlessly, at the touch of a button.

Intuitive assembly

Easy visual matching of parts for intuitive assembly.

Fully compatible range

This breast pump can be used in combination with other feeding products in the Philips Avent range, including our Classic bottles and milk storage containers. Philips Avent also provides a range of breast care accessories to help you breastfeed for longer and enhance your comfort.

Specifications

Breast pump

Material: BPA-free*

Breastfeeding support service

Professional support: unlimited one-to-one video calls with a breastfeeding professional, advice for preparing for breastfeeding and what to expect when breastfeeding a newborn, advice on breastfeeding issues

Design

Bottle design: Ergonomic shape, Wide neck

Breast pump design: Compact design

Country of origin

England

Material

Bottle: Polypropylene, BPA-free*

Breast pump: Polypropylene, BPA-free*

Teat: Silicone, BPA-free*

What is included

Breast pad sample packs: 2 (2 Day pads and 2 Night pads) pcs

Extra-soft newborn flow teat: 1 pcs

Sealing disc for milk storage: 1 pcs

Natural bottle 125 ml/4 oz: 1 pcs

Breast pump body: 1 pcs

Standard-size cushion: 1 pcs

Base unit incl. tubing: 1 pcs

Spare diaphragm: 1 pcs

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 3 Expression settings, 1 Stimulation mode

Development stages

Stage: 0–6 months

