



Philips Avent  
Single electric breast pump

**Ultra Comfort**  
Includes 4oz Natural bottle



SCF332/61

## More comfort, more milk\*

Clinically proven comfort\*

Sit comfortably without leaning forward as our soft massage cushion gently stimulates your milk flow. Use our quiet pump anytime, anywhere, even with batteries. It's easy to set up, personalize, use and clean.

### Effective expression

- 3 customizable settings at your fingertips

### Easy to set up, express, clean and travel with

- Easy to use anywhere
- Intuitive design, easy to assemble and clean
- No spills or messy transfers

### Express in comfort

- Clinically proven less leaning forward\*
- Designed to let you express quietly
- Soft cushion with massaging petals

### The support in your breastfeeding journey

- Compatible with our natural bottles, for easy feeding
- Easy access to warranty and support

**PHILIPS**  
**AVENT**

# Highlights

## Relaxing expression position



Research\* suggests that when you're relaxed, your milk flows more easily. Thanks to its unique design, you won't need to learn forward when you pump. With this breast pump, you can comfortably sit up while your milk flows directly from your breast into the bottle.

## Soft massage cushion



The massage cushion is made of a soft, velvety material to give you a feeling of warmth and helps stimulate your milk flow. With nature in mind, the iconic petal cushion is designed to mimic baby suckling, which is intended to gently help stimulate let down.

## Quiet pump

Designed for discreet expression, anytime.

## Control breast milk flow

With the press of a button, you can start pumping in the gentle stimulation mode to activate your let-down. To begin breast milk collection, choose one of the three flow settings for your preferred suction speed. To

increase the suction, press a button again to adjust.

## Compact and portable

Our expression kit, which includes the funnel, cushion and bottle, is easy to hold in one hand and position on your breast. Place the lightweight base unit within easy reach for full control when pumping. It's also very convenient when you go out. Just wrap the tube around the base unit. If you plan to travel abroad, read the specifications tab to check the voltage before plugging it in or use the batteries for easy back-up.

## Intuitive assembly



The breast pump features an intuitive design that is easy to clean. And since it only has a few parts, you can quickly assemble it and get started. Before each use, we recommend you wash the parts in warm, soapy water and rinse. Then sanitize with a Philips Avent steam sterilizer, or boil them for five minutes for best results. Use a damp cloth to clean the motor unit, cap, silicone tube and adapter.

## No milk spillage in the tubes



Thanks to its closed system, the breast pump collects your milk directly into the storage container. This way your milk is collected while the suction tubes remain clean.

## Handy accessories



To help you and baby on your breastfeeding journey, you'll find a 4oz Philips Natural bottle included with your pump. The wide, breast-shaped teat promotes natural latch-on similar to the breast to help make it easier for your baby to combine breast with bottle feeding. The pump also fits our breast feeding accessories, such as milk storage cups. See our Accessories tab for more options and replacement parts.

## Warranty coverage

Your breast feeding journey is an ongoing experience. That's why we offer specialized product warranty coverage. Register your breast pump for product coverage as well as special offers and easy access to product support. Click the Support tab to see eligible details in your country.

# Specifications

## Design

- Breast pump design: Compact design
- Bottle design: Ergonomic shape, Wide neck

## Country of origin

- England

## Material

- Bottle: BPA free\*, Polypropylene
- Nipple: BPA free\*, Silicone
- Breast pump: BPA free\* (food contact parts only)

## What is included

- Sealing disc for milk storage: 1 pcs
- Natural bottle 4oz: 1 pcs
- Breast pump body: 1 pcs
- Base unit incl. tubing: 1 pcs

- Travel cover: 1 pcs
- Standard size cushion(19.5 mm): 1 pcs (a larger size cushion is available and sold separately)
- Disposable breast pads: 4 pcs

## Ease of use

- Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

## Functions

- No leaning forward: Sit in a comfortable position
- Soft massage cushion: Gentle stimulation
- Settings: 1 Stimulation mode, 3 Expression settings

## Development stages

- Stage: 0 - 6 months



Issue date 2023-10-05

Version: 3.3.0

© 2023 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* Clinically proven comfort: In a test carried out on 110 mothers in the USA, UK, China and Russia in March 2016, mothers awarded an average of 8.6/10 score for Philips Avent performance on comfort.

\* More milk: Independent research has shown that there may be a link between higher comfort levels and increased milk production. See [www.philips.com/AVENT](http://www.philips.com/AVENT)

\* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011