

PHILIPS

Handle for breast
pump

ISIS

manual breast pump

SCF923



To make the breast pump diaphragm move

Find matching products on the specifications tab

Pushing this handle, a part of your Avent manual breast pump, up and down creates suction. You can control the suction by varying the pressure on the handle.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCF310/08, SCF310/60, SCF300/13, SCF300/20, SCF290/12, SCF290/20, SCF300/12, SCF310/20, SCF290/13, SCF310/13, SCF310/12, SCF300/60, SCF300/08

