

PHILIPS

12mm Micro-
Flexible Tubing

DreamStation Go

Philips smallest tubing

Offers freedom of movement

Philips lightest tubing



HH1455



Freedom of movement when you sleep

Our 12mm tubing is 58% smaller than Philips 22mm performance tubing* and offers better flexibility and ease of movement when compared to a 15mm tubing.**

The 12mm tubing is our smallest and lightest ever

- Smallest and lightest

Highlights

Smallest and lightest

The 12mm micro-flexible tubing is the smallest and lightest tube Philips has ever designed. In a recent head-to-head trial, patients found the DreamStation 12mm micro-flexible tubing to

offer better flexibility, more comfort, and greater ease of movement when compared to their standard tubing.**

Specifications

Technical information

Material: TPE and Polypropylene

Replacement: Every 6 months

Tubing length: 6.15ft, 187.5 cm

Tubing weight: 2.44 Oz, 69 grams

Tubing inner diameter: 12.12mm

Part Number

Micro-flexible 12mm tubing: PR12

Service

Warranty: 90-days

Compatibility

PAP System: DreamStation Go CPAP, AutoPAP



* **When compared to Philips Respironics 22mm performance tubing

** **Head-to-head, incentivized trial, of DreamStation 15mm tubing vs. DreamStation 12mm micro-flexible tubing, 2016 (n=31)