

# PHILIPS

## sonicare



Sonic electric  
toothbrush

4100 Series

Sonicare technology  
QuadPacer and SmarTimer  
Slim ergonomic design  
14-day battery life



HX3689/24

## Goodbye manual toothbrush. Hello Sonicare.

Up to 7X better plaque removal\*

Sonicare technology combined with our brush action gently removes up to 7x better\* plaque than a manual toothbrush.

### **Designed around you**

- Use less. Reuse more
- Always know when to replace your brush heads
- 14-day long lasting battery life
- Switch between 2 intensity settings

### **Optimize your brushing**

- Optimize your brushing with SmarTimer and QuadPacer

### **Protects your gums**

- Pressure Sensor helps protect your teeth and gums

### **Provides superior clean**

- Our unique technology brings you a powerful yet gentle clean
- Up to 7x better plaque removal vs a manual toothbrush\*

# 4100 Series

Sonic electric toothbrush

HX3689/24

## Specifications

### Power

Voltage: DC5V

### Technical specifications

Operating time (full to empty): 14 days\*\*\*

Battery: Rechargeable

Battery type: Lithium ION

### Design and finishing

Color: Lush Lilac Prismatic

### Service

Warranty: 2-year limited warranty

### Ease of use

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

**Brush head system:** Easy click-on brush heads

**BrushSync Replacement Reminder:** Replacement reminder, Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

### Items included

**Handle:** 1 4100 Series

**Brush heads:** 1 C2 Optimal Plaque Control

**Charger:** 1

### Cleaning performance

**Performance:** Removes up to 7x more plaque\*

**Timer:** Quadpacer and SmarTimer

**Pressure feedback:** Vibrates handle to alert user

### Modes

**Clean:** For exceptional everyday clean

**2 intensities:** High, Low

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**BrushSync Replacement Reminder:** Always know when to, replace brush heads

