

# IHE Integration Statement forConnect



| Vendor  | Product Name   | Version                                     | Date       |
|---|--|---|------------|
| Forcare   | forConnect   | 2023-2                                      | 2023-10-12 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below |  |   |            |
| Integration Profile   | Actor  | Option                                      |            |
| Consistent Time (CT)  | Time Client (ITI-1)                                  |   |            |
| Audit Trail and Node Authentication (ATNA)  | Secure Application (ITI-19, ITI-1, ITI-20)           |   |            |
| Cross-Enterprise User Assertion (XUA)   | X-Service-User (ITI-40)                              |   |            |
| Internet User Authorization (IUA)   | Resource Server (ITI-72)                             | SAML Token                                  |            |
| Cross-Enterprise Document Sharing (XDS.b)   | Document Source (ITI-41)                             | Document Addendum                           |            |
|   |  | Document Replacement                        |            |
|   | Document Consumer (ITI-18, ITI-43)                   | Document Transformation                     |            |
|   |  | Document Metadata Update                    |            |
| Document Administrator  | On-Demand Documents                                  |   |            |
|   | Delayed Document Assembly                            |   |            |
| Cross-Enterprise Document Reliable Interchange (XDR)  | Document Source (ITI-41)                             | Update Metadata (ITI-57)                    |            |
| Cross-Enterprise Document Sharing for Imaging (XDS-I)   | Imaging Document Source (ITI-15)                     | Transmit Home Community ID                  |            |
|   |  | Set of DICOM Instances                      |            |
|   |  | PDF Report                                  |            |
| Imaging Document Consumer (RAD-55)  | Text Report  |   |            |
|   | Set of DICOM Instances (RAD-55, RAD-69)              |   |            |
|   | PDF Report   |   |            |
| Cross-Enterprise Document Sharing for Imaging (XDS-I.b)   | Imaging Document Source (ITI-41, RAD-68)             | CDA Wrapped Text Report                     |            |
|   |  | CDA Imaging Report with Structured Headings |            |
|   |  | Imaging Document Consumer (RAD-55, RAD-69)  |            |
| Cross Community Access (XCA)  | Responding Gateway (ITI-38, ITI-39)                  |   |            |
| Patient Identifier Cross-referencing for MPI (PIX)  | Patient Identifier Cross-Reference Consumer (ITI-9)  |   |            |
| Patient Identifier Cross-referencing HL7V3 (PIXV3)  | Patient Identifier Cross-Reference Consumer (ITI-45) |   |            |
| Basic Patient Privacy Consents (BPPC)   | Content Creator                                      |   |            |
| Scanned Documents (XDS-SD)  | Content Creator                                      |   |            |
| Cross-Enterprise Document Workflow Profile (XDW)  | Content Creator                                      | XDW Cross-Community with XDR                |            |
|   | Content Consumer                                     | Document Import Option                      |            |
|   | Content Updater                                      | Document Import Option                      |            |
|   |  | XDW Cross-Community with XDR                |            |
| Cross Enterprise Basic eReferral Workflow Definition Profile (XBeR-WD)  | Referral Requester                                   |   |            |
| Internet address for vendor's IHE information <a href="http://www.forcare.com">www.forcare.com</a>  |  |   |            |

## Links to Standards Conformance Statements for the Implementation

HL7 [www.forcare.com](http://www.forcare.com)

---

DICOM [www.forcare.com](http://www.forcare.com)

---

**Links to general information on IHE**

---

In North America: [www.ihe.net](http://www.ihe.net)

In Europe: [www.ihe-europe.net](http://www.ihe-europe.net)

In Japan: [www.ihe-j.org](http://www.ihe-j.org)

---

© 2006-2023 Forcare Holding B.V. All rights reserved.

---