

# IHE Integration Statement forNotify



Vendor	Product Name	Version	Date
Forcare	forNotify	2023-2	2023-10-12

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below

Integration Profile	Actor	Option
Consistent Time (CT)	Time Client	
Audit Trail and Node Authentication (ATNA)	Secure Application (ITI-19, ITI-1, ITI-20)	
Cross-Enterprise User Assertion (XUA)	X-Service-Provider (ITI-40)	Subject-Role
		PurposeOfUse
	X-Service-User (ITI-40)	Subject-Role
		PurposeOfUse
Internet User Authorization (IUA)	Resource Server (ITI-72)	SAML Token
Document Metadata Subscription (DSUB)	Document Metadata Notification Recipient (ITI-53)	
	Document Metadata Notification Subscriber (ITI-52)	
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer (ITI-18, ITI-43)	
Cross-Enterprise Document Workflow (XDW)	Content Consumer	Document Import Option

Internet address for vendor's IHE information [www.forcare.com](http://www.forcare.com)

## Links to Standards Conformance Statements for the Implementation

HL7	<a href="http://www.forcare.com">www.forcare.com</a>
DICOM	<a href="http://www.forcare.com">www.forcare.com</a>

## Links to general information on IHE

In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>
--	---	--

© 2006-2023 Forcare Holding B.V. All rights reserved.