



Philips Minicells  
Battery

1.4V  
Zinc-air

ZA13BX30



# Top quality Zinc-air technology for hearing aids

Provides a long lasting performance in personal hearing aids

## Top Performance

- Reliable for long term performance
- Zinc-air technology is designed for hearing aids
- The battery remains fresh for use for up to 3 years

## Ease of Use

- Colors make recognition of battery sizes easier
- The dispenser makes removal of cell easy

# PHILIPS

Battery  
1.4V Zinc-air

# Specifications

## Power

- Battery type: Hearing Aid / Zinc-Air
- Battery voltage: 1.4 V

## Green Specifications

- Heavy metals: Cd free, Pb free
- Packaging material: Carton, PET
- Packaging type: PET blister

## Technical specifications

- Shelf life: 3 years

## Outer Carton

- Outer carton (L x W x H): 28 x 21.5 x 18 cm
- Outer carton (L x W x H): 11.0 x 8.5 x 7.1 inch
- Nett weight: 2.2 kg
- Nett weight: 4.850 lb
- Gross weight: 3.25 kg
- Gross weight: 7.165 lb
- Tare weight: 1.05 kg

- Tare weight: 2.315 lb
- GTIN: 1 48 95229 12224 3
- Number of consumer packagings: 50

## Packaging dimensions

- Packaging dimensions (W x H x D):  
5 x 9.6 x 2.8 cm
- Packaging dimensions (W x H x D):  
2.0 x 3.8 x 1.1 inch
- Nett weight: 0.044 kg
- Nett weight: 0.097 lb
- Gross weight: 0.055 kg
- Gross weight: 0.121 lb
- Tare weight: 0.011 kg
- Tare weight: 0.024 lb
- EAN: 48 95229 12224 6
- Type of shelf placement: Both
- Number of products included: 6
- Packaging type: Blister

ZA13BX30/00

# Highlights

## Reliable performance

In manufacturing excellent technology is being used to guarantee a reliable and long term performance.

## Zinc-air technology

Zinc-air technology to power your hearing aid.

## Easy color coding

To make recognition of the most common battery sizes easier, we use color coding for battery sizes according to industry standards.

## Easy rotary dispenser

In this packaging system a rotary wheel is used which makes removal of the cell very easy.



Issue date 2022-06-22

Version: 2.1.2

12 NC: 8670 001 77063  
EAN: 48 95229 12224 6

© 2022 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)