



Battery

- 1.4V
- Zinc-air

ZA13B6A/59

Top quality Zinc-air technology for hearing aids

Provides a long lasting performance in personal hearing aids

Benefits

Top Performance

- Reliable for long term performance
- Zinc-air technology is designed for hearing aids
- The battery remains fresh for use for up to 3 years

Ease of Use

- Colors make recognition of battery sizes easier
- The dispenser makes removal of cell easy

Features

Reliable performance

In manufacturing excellent technology is being used to guarantee a reliable and long term performance.

Zinc-air technology

Zinc-air technology to power your hearing aid.

Easy color coding

To make recognition of the most common battery sizes easier, we use color coding for battery sizes according to industry standards.

Easy rotary dispenser

In this packaging system a rotary wheel is used which makes removal of the cell very easy.

Specifications

[field not found 'disclaimer_text']

Power

Battery type

Hearing Aid / Zinc-Air

Battery voltage

1.4 V

PHILIPS

Minicells

Green Specifications

Heavy metals	Cd free Pb free
Packaging material	Carton
Packaging type	PET PET blister

Technical specifications

Shelf life	3 years
------------	---------

Outer Carton

GTIN	1 48 95229 12558 9
Number of consumer packagings	200
Gross weight	2.1 kg
Height	13 cm
Length	23 cm
Nett weight	1.00800 kg

Tare weight	1.09200 kg
Width	22 cm

Packaging dimensions

EAN	48 95229 12558 2
Number of products included	6
Depth	0.7 cm
Gross weight	0.0086 kg
Height	9.4 cm
Nett weight	0.00504 kg
Tare weight	0.00356 kg
Width	4.77 cm

Product dimensions

Depth	0.79 cm
Height	0.54 cm
Width	0.79 cm



data subject to change
2023, June 21

Version: 1.1.2
EAN: 4895229125582

© 2023 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com