

# PHILIPS

## Battery

## Minicells

10/230

Zinc-air

ZA10B6A



## Top quality Zinc-air technology for hearing aids

Provides a long lasting performance in personal hearing aids

### **Ease of Use**

- Colors make recognition of battery sizes easier
- The dispenser makes removal of cell easy

### **Top Performance**

- Reliable for long term performance
- Zinc-air technology is designed for hearing aids
- The battery remains fresh for use for up to 3 years

## Highlights

### Reliable performance

In manufacturing excellent technology is being used to guarantee a reliable and long term performance.

### Zinc-air technology

Zinc-air technology to power your hearing aid.

### Easy color coding

To make recognition of the most common battery sizes easier, we use color coding for battery sizes according to industry standards.

### Easy rotary dispenser

In this packaging system a rotary wheel is used which makes removal of the cell very easy.

## Specifications

### Green Specifications

**Heavy metals:** Cd free, Pb free

**Packaging material:** Carton, PET

**Packaging type:** PET blister

### Outer Carton

**Length:** 22 cm

**Number of consumer packagings:** 200

**Width:** 21 cm

**Gross weight:** 1.4 kg

**Height:** 13 cm

**GTIN:** 1 48 95229 12557 2

**Nett weight:** 0.33600 kg

**Tare weight:** 1.06400 kg

### Power

**Battery Capacity:** max. 90 mAH

**Battery type:** Hearing Aid / Zinc-Air

**Battery voltage:** 1.4 V

### Packaging dimensions

**Height:** 9.4 cm

**Width:** 4.77 cm

**Depth:** 0.5 cm

**Number of products included:** 6

**EAN:** 48 95229 12557 5

**Gross weight:** 0.0054 kg

**Nett weight:** 0.00168 kg

**Tare weight:** 0.00372 kg

### Dimensions

**Master carton quantity:** 100

**Product dimensions (WxDxH):** 0.36 x 0.58 cm

**Outer carton dimensions WxDxH:** 10.2 x 17 x 11.5 cm

**Inner carton dimensions WxDxH:** 3.2 x 5 x 9.8 cm

**Product weight:** 0.0002 kg

### Technical specifications

**Shelf life:** 3 years

