



Philips Minicells
Battery

Zinc-air

ZA13BX60



Top-quality Zinc-air technology for hearing aids

Provides long-lasting performance in personal hearing aids

Top Performance

- Reliable for long-term performance
- Zinc-air technology is designed for hearing aids
- The battery remains fresh for use for up to 3 years

Ease of Use

- Colours make recognition of battery sizes easier
- The dispenser makes removal of cell easy

PHILIPS

Battery
Zinc-air

ZA13BX60/00

Specifications

Power

- Battery type: Hearing Aid/Zinc-Air
- Battery voltage: 1.4 V

Green Specifications

- Heavy metals: Cd free, Pb free
- Packaging material: Carton, PET
- Packaging type: PET blister

Technical specifications

- Shelf life: 3 years

Outer Carton

- GTIN: 1 48 95229 12225 0
- Number of consumer packages: 60

Packaging dimensions

- EAN: 48 95229 12225 3
- Number of products included: 6

Highlights

Reliable performance

Excellent technology is used in manufacturing to guarantee reliable and long-term performance.

Zinc-air technology

Zinc-air technology to power your hearing aid.

Easy colour-coding

To make it easier to recognise the most common battery sizes, we use colour-coding for battery sizes according to industry standards.

Easy rotary dispenser

A rotary wheel is used in this packaging system, which makes removal of the cell very easy.



Issue date 2022-05-04

Version: 1.1.1

12 NC: 8670 001 77064
EAN: 48 95229 12225 3

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com