




Philips Minicells
Battery

10/230
Zinc-air

ZA10B6A



Top quality Zinc-air technology for hearing aids

Provides a long lasting performance in personal hearing aids

Top Performance

- Reliable for long term performance
- Zinc-air technology is designed for hearing aids
- The battery remains fresh for use for up to 3 years

Ease of Use

- Colors make recognition of battery sizes easier
- The dispenser makes removal of cell easy

PHILIPS

Specifications

Power

- Battery type: Hearing Aid / Zinc-Air
- Battery voltage: 1.4 V
- Battery Capacity: max. 90 mAh

Green Specifications

- Heavy metals: Cd free, Pb free
- Packaging material: Carton, PET
- Packaging type: PET blister

Technical specifications

- Shelf life: 3 years

Dimensions

- Inner carton dimensions WxDxH: 3.2 x 5 x 9.8 cm
- Master carton quantity: 100
- Outer carton dimensions WxDxH: 10.2 x 17 x 11.5 cm

- Product dimensions (WxDxH): 0.36 x 0.58 cm
- Product weight: 0.0002 kg

Outer Carton

- GTIN: 1 48 95229 12557 2
- Number of consumer packagings: 200
- Gross weight: 1.4 kg
- Outer carton (L x W x H): 22 x 21 x 13 cm
- Nett weight: 0.33600 kg
- Tare weight: 1.06400 kg

Packaging dimensions

- EAN: 48 95229 12557 5
- Number of products included: 6
- Packaging dimensions (W x H x D): 4.77 x 9.4 x 0.5 cm
- Gross weight: 0.0054 kg
- Nett weight: 0.00168 kg
- Tare weight: 0.00372 kg

Highlights

Reliable performance

In manufacturing excellent technology is being used to guarantee a reliable and long term performance.

Zinc-air technology

Zinc-air technology to power your hearing aid.

Easy color coding

To make recognition of the most common battery sizes easier, we use color coding for battery sizes according to industry standards.

Easy rotary dispenser

In this packaging system a rotary wheel is used which makes removal of the cell very easy.



Issue date 2022-05-12

Version: 1.1.2

12 NC: 8670 001 80734
EAN: 48 95229 12557 5

© 2022 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com