

PHILIPS

sonicare

Standard sonic
toothbrush heads

G3 Premium Gum
Care

4-pack

Standard size

Click-on

BrushSync mode pairing

HX9054/17



Our brush head for superior gum health

Up to 7 x healthier gums in just 2 weeks*

Healthy teeth start with healthy gums. Take care of yours with our G3 Premium Gum Care brush head. Soft sides flex and adapt to the contours of your gums, giving you 2 x more surface contact** for superior plaque removal at the gum line.

Innovative technology

- Philips Sonicare's advanced sonic technology
- Automatically selects optimal mode for ultimate results****

Designed to optimise performance

- Always know when to replace. Always get an effective clean.
- Click-on design for simple brush head placement
- Tested to meet your oral health needs

For exceptional gum care

- Up to 7 x healthier gums in just 2 weeks***

An unrivalled deep clean

- Soft rubber sides flex and provide up to 2x more surface contact

Now more sustainable

- 75% of the plastic in this product is bio-based*****

G3 Premium Gum Care

Standard sonic toothbrush heads

HX9054/17

Specifications

Design and finishing

Colour: White

Bristle stiffness feel: Soft

Reminder bristles: Blue bristle colour fades away

Material brush head: Soft, flexible rubber sides

Size: Standard

Smart Brush Head Recognition

Compatibility

Brush head system: Click-on

Suitable for these models: 2 Series plaque defence, 2 Series plaque defence, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp

Not suitable for: Philips One

Items included

Brush heads: 4 G3 Premium Gum Care

Quality and performance

Replacement: Every 3 months

Tested: for optimal usage

Health benefits

Gum health: Up to 7 x healthier gums*

Plaque removal: Removes 10 x more plaque*****

