

# PHILIPS



Sonic electric  
toothbrush with app

DiamondClean  
Smart 9750

Up to 20x more plaque removal\*

Up to 15x healthier gums\*\*

Up to 100% less stains\*\*\*

Visible pressure alert



HX9924/67

## Complete care for a healthier mouth

Feel confident with 100% coverage

The DiamondClean Smart provides a deep clean for a whiter, healthier smile. High-performance brush heads let you focus on all areas of your oral health and our Smart Sensor technology gives you personalized feedback and coaching.

### Proven to improve oral health

- Enjoy long-lasting fresh breath

### Personalized in-app coaching to perfect brushing

- Smart guidance for gentle care
- Guides you and provides instant feedback as you brush
- Helps you improve your brushing technique
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas

### Optimize your brushing

- Make sure you're getting the most from your brush head
- Automatically selects the best mode for you
- Choose from 5 modes and 3 intensity settings
- Easy and stylish charging - home and away

### Premium performance

- Up to 7x healthier gums in just 2 weeks\*\*\*\*\*

# DiamondClean Smart 9750

Sonic electric toothbrush with app

HX9924/67


## Highlights

### 5 modes, 3 intensity settings




DiamondClean Smart comes with five modes which cover a range of brushing needs: Clean mode is for exceptional daily cleaning, White+ is for removing stains, Deep Clean+ gives you an invigorating deep clean, Gum Health provides a gentle yet effective clean for gums, and Tongue Care is the perfect mode for cleaning your tongue. It breaks up tongue debris and sweeps away bad-breath bacteria. Three intensities allow you to choose between a higher setting to boost your clean and a lower one for more-sensitive mouths.

### TouchUp feature




If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

### Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

### TongueCare+ tongue brush



Click on the TongueCare+ tongue brush to gently remove odor causing bacteria from the pores of your tongue. Its 240 specially designed MicroBristles get between all your tongue's ridges and grooves to remove bacteria and debris that lead to bad breath. Team with our antibacterial BreathRx tongue spray for superior cleaning and super fresh breath.

# DiamondClean Smart 9750

Sonic electric toothbrush with app

HX9924/67

## Specifications

### Connectivity

**Bluetooth® wireless technology:** Connected brushing app

### Power

**Voltage:** 110-220 V

### Technical specifications

**Battery:** Rechargeable

**Operating time (full to empty):** 14 days\*\*\*\*\*

**Battery type:** Lithium ION

### Design and finishing

**Color:** Rose gold

### Service

**Warranty:** 2-year limited warranty

### Compatibility

**Android compatibility:** Android phones with Bluetooth 4.0 and above, Search Sonicare on Google Play Store for detailed OS compatibility.

**iOS compatibility:** iPhone with Bluetooth 4.0 and above, Search Sonicare on App Store for detailed iOS compatibility.

### Ease of use

**Display:** Illuminated display

**Battery indicator:** Illuminated icon indicates battery life

**Handle:** Slim ergonomic design

**Handle compatibility:** Easy click-on brush heads

**Brush head recognition:** Syncs with optimal mode

**Replacement reminder:** To always ensure best results, reminder icon lights up

### Items included

**Handle:** 1 DiamondClean Smart

**Travel case:** 1 USB charging case

**Brush heads:** 1 A3 Premium All-in-One, 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 TongueCare+ tongue brush

**Glass charger:** 1

### Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Health benefits:** Up to 15x healthier gums\*\*

**Performance:** Removes up to 20x more plaque\*

**Timer:** BrushPacer and SmarTimer

**Whitening benefits:** Up to 100% less stains\*\*\*

**Pressure feedback:** Ring lights up in purple, Vibrates handle to alert user

**TouchUp:** Ensures 100% coverage

### Modes

**3 intensities:** High, Medium, Low

**Clean:** For exceptional everyday clean

**Deep Clean+:** For an invigorating deep clean

**Gum Health:** Special attention to molars

**TongueCare+:** For long lasting clean breath

**White+:** To remove surface stains

### BrushSync Mode Pairing

**C3 Premium Plaque Control:** Pairs with Clean mode

**G3 Premium Gum Care brush head:** Pairs with Gum Health mode

**TongueCare+ tongue brush:** Pairs with TongueCare mode

**A3 Premium All-in-One:** Pairs with Clean mode

### Smart sensor technology

**Pressure sensor:** Alerts when brushing too hard

**Location sensor:** Tracks & improves coverage

**Scrubbing sensor:** Guides to reduce scrubbing

**3D mouth map:** Highlights trouble areas

**TouchUp:** To address missed spots

**Position sensor:** Guides to areas you missed



\* vs. a manual toothbrush.

\*\* \*when used with A3 brush head vs. a manual toothbrush.

\*\*\*vs. a manual toothbrush. In a lab study, actual results may vary.

\*\*\*\*than a DiamondClean brush head.

\*\*\*\*\*in Gum Health mode vs a manual toothbrush; measured by GBI

\* \*\*\*\*\*based on two periods of two-minute brushings per day on standard mode.