

**PHILIPS**

**AVENT**



## Giftset

### Manual Breast Pump

All-in-one set

Sterilise, pump and feed



SCD430/50

## Sterilise, pump and feed

### Manual breast pump set with steriliser

This Breastfeeding set includes everything you need from pumping, feeding and cleaning up. It includes a variety of accessories: a manual breast pump, steriliser, soothers and bottles. Finding the right teat is important. See more below

#### **Other benefits**

- Choose the right teat flow for your baby

#### **Feed with Philips Avent bottles**

- For a natural latch-on

#### **Sooth with ultra soft soothers**

- Flexible shield for a comfortable fit

#### **Express with the Manual Breast Pump**

- For a quick milk flow
- More comfortable experience
- Ideal for occasional pumping any time, anywhere
- No need to lean forwards

#### **Sterilise with the Microwave steriliser**

- Fast and effective
- Contents stay sterile for up to 24 hours

## Highlights

### Natural Motion Technology



Enjoy the gentle comfort of the Philips Avent portable manual breast pump. With Natural Motion Technology, inspired by baby's own suckling action for a quick letdown. You can easily adjust the rhythm and vacuum.

### Adaptive silicone cushion



The silicone cushion gently flexes and adapts to fit your unique nipple's size and shape. Fits 99.98% of all nipple sizes.

### Finding the right teat is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a teat with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

### Find the right teat flow




Every baby feeds differently and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalise your bottle. All of the Natural Response Teats are made of soft silicone.

### Portable and easy to use



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on the go simple and discreet.

### Sit comfortably



Relax and sit comfortably while you express - thanks to our design that lets you sit upright, rather than having to lean forward.

### Sterilise in just 2 minutes



The microwave steam steriliser can be used to sterilise baby bottles and other products in just 2 minutes, killing 99.9% of germs and bacteria. The exact length of the cycle depends on the wattage of your microwave, of course. 2 minutes at 1200-1850 W, 4 minutes at 850-1100 W, 6 minutes at 500-800 W

# Manual Breast Pump

Giftset

SCD430/50

## Highlights

**Kills 99.9% of germs**



The microwave steam steriliser is proven to kill 99.9% of germs and bacteria.

**Includes our Natural bottles**



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

**Includes ultra soft soothers**



Take care of sensitive skin with the Philips Avent ultra soft soother. Our flexible shield follows the curves of your baby's cheeks, leaving fewer marks and less irritation for more comfortable soothing.

# Manual Breast Pump

Giftset

SCD430/50

## Specifications

### What is included

Bottle brush: 1 pcs  
Manual breast pump: 1 pcs  
Microwave steriliser: 1 pc  
ultra soft white 0-6 m soothers: 2 pcs

125 ml (4 oz) Bottle: 2 pcs  
260 ml (9 oz) bottles: 2 pcs  
Slow Teat: 2 pcs  
Medium Teat: 2 pcs

### Functions

**Soft and adaptive cushion:** Gentle stimulation  
**Portable and lightweight:** Easy on the go

© 2025 Koninklijke Philips N.V.  
All Rights reserved.

Issue date 2025-04-15  
Version: 8.8.1

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

EAN: 87 20689 00811 2

[www.philips.com](http://www.philips.com)



\* Based on: (1)Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index and breast anatomy: are they related?. Breastfeeding Medicine, 26 April 2019, (109 participants, Israel); (2)Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

\* Journal of Obstetric, Gynecologic & Neonatal Nursing, May 1993, (20 participants (Caucasian), USA); (3)Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia)

\* BPA-free breast pump: Only associated with the bottle and other parts that come into contact with breast milk. Following EU regulation, 10/2011