

PHILIPS

Mask - Fit Pack (S,
M, L)

DreamWear

3 cushion sizes

Medium frame



HH1124



Like wearing nothing at all*

With the freedom and comfort of DreamWear Gel Pillows, you can sleep how you want.

Open field of vision

- Open design offers an extremely wide field of vision

Innovative design

- Air tube connection at the top of the head
- Gel pillows cushion
- Modular design
- Supports your sleep therapy

Comfortable design

- Easy to fit
- Slim-line headgear
- Soft fabric wraps
- Soft, flexible frame

Highlights

Easy to fit



Choose the cushion so the pillows cushion blue gel bases sit comfortably under the nose without any gaps. Only the soft silicone tip should be inside your nostrils. Do not push any part of the blue gel base inside your nose.

Sleep how you want



By moving the tube location to the top of the head, DreamWear Gel Pillows offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear Gel Pillows.

Gel pillows cushion




Less discomfort on the face, nose and nostrils. Users rated DreamWear Gel Pillows more comfortable than their prescribed mask.

Modular design



Modular design for quick assembly and cleaning, with interchangeable pillows and nasal cushions.

Slim-line headgear



DreamWear combines comfort with aesthetics. Users reported that the DreamWear Gel Pillows mask came closest to making them feel like they have nothing on their face versus their prescribed mask.

Soft fabric wraps




Provides a soft interface to the cheeks.

Soft, flexible frame




The flexible material through which the air flows is applied gently to the skin and is pleasant to wear.

Unobstructed view



Open design makes it easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Supports your sleep therapy



DreamWear gives you the freedom to choose the cushion that works best for you, giving you freedom of movement and better quality of sleep.

Specifications

Package contents

Include: Small Cushion, Medium Cushion, Large Cushion, Medium Frame

Intent of use

To be used on patients <math>< /math>< /math>>30 kg: with CPAP prescription

Cleaning: Headgear, tubing

Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Cleaning: Mask

Hand wash with mild detergent: Rinse thoroughly, air dry¹

Replacement

Inspect daily for wear: Replace at first sign of wear

Product details

Materials: Nylon, Polyurethane foam, Polycarbonate, Polyester, Polyurethane gel, Silicone, Spandex

Warranty: 90 days

Operating pressure: 4 to 20 cmH₂O

Sound level: 18 dBA

Number of parts: 5



* 2017 Philips User Preference Questionnaire

* The scope of supply includes a standard mask frame (size M) and three cushions for all sizes (S, M, L). The sleep therapy device is not included.

* The costs of DreamWear are private and are usually not reimbursed by the statutory health insurance.