

# PHILIPS

Compressor  
Nebulizer System

InnoSpire Essence

Compact home-use device

Fast, effective treatment

Nebulizer kit included



Product Picture

HH1340



## A nebulizer system you can rely upon

When reliability is key choose InnoSpire Essence. Essence provides fast, efficient nebulized drug therapy for you and your family at home.

### Fast treatment times

- Boosts air flow for a faster treatment time

### Proven SideStream technology

- Proven nebulizer technology
- Robust design meets reliability

### Short treatment times

- Treatment delivered in 6-8 minutes\*

### For use in the home

- Designed for intermittent use
- Compatible with commonly prescribed nebulizer medications

# Highlights

## For home use

Designed for intermittent use and suitable for the home, InnoSpire Essence provides the consistent delivery high-quality aerosol treatment that you can rely on.

## Treatment time

InnoSpire Elegance combined with the SideStream Nebulizer, is designed to provide fast drug delivery and short treatment times.

## Versatility

Capable of nebulizing a wide variety of medications that are routinely with nebulizers used to treat various respiratory conditions.

## SideStream nebulizer system



Together with the SideStream reusable nebulizer, InnoSpire Elegance provides consistently fast and efficient aerosol therapy.

## Proven nebulizer technology

Reduced Treatment Times - Our Sidestream nebulizer is specially designed to keep your treatment times as short as possible. As the InnoSpire Essence pushes air through the tube, additional air is drawn in through the top of the nebulizer to create a fast flow of aerosol medication for you to inhale. Treatments can take as little as 6 minutes\* Consistent delivery of high quality aerosol therapy - Sidestream is designed to withstand frequent and repeated use and will continue to efficiently deliver your nebulizer medication until it is time for replacement.

## Power and style all in one

InnoSpire Essence is designed to both fit your lifestyle and cope with frequent use. You can be confident that you will receive high quality aerosol for the life of the nebulizer.

# Specifications

## Design

**Compact:** Weighs only 1.5 kg

## Operating parameters

**Electrical rating/power:** 230 V / 50 Hz, 1 Amp

**Maximum flow rate:** 8 LPM

**Mode of operation:** Intermittent (30 min. on/off)

**Average flow rate:** 6 LPM @ 10 psi

## Performance data

**MMAD:** 2.9 µm

**Fine particle fraction:** 77% below 5 µm

**Nebulisation capacity:** 8 ml

## Service

**Warranty:** 3 years

## Product details

**Weight:** 1.5 kg

**Size:** 165x165x108 mm

**Noise level:** 58 ±3 dBA

**Treatment time (6-8 minutes\*):** \* 2.5ml Salbutamol dose

