

# PHILIPS

## AVENT

Natural Response  
Baby Bottle

4 Bottles

9oz/260ml

Flow 3 nipple

1m+



SCY903/68



## Supports baby's individual drinking rhythm A nipple that works like a breast

The Natural Response Nipple releases milk only when baby actively drinks. Babies can drink at their own rhythm like on the breast, so easy to combine breast and bottle feeding. Finding the right nipple is important. See more info below.

### Easy to combine breast and bottle feeding

- Natural latch on with breast-shaped nipple
- Nipple releases milk when baby actively drinks

### Other benefits

- Choose the right nipple flow for your baby
- No-drip nipple design prevents spills and lost milk
- Compatible across the Philips Avent range
- Easy to hold even for little hands
- Simple to use, easy to clean and quick to assemble
- Natural Response Nipples and Bottles are BPA free\*

### Vents air away from your baby's tummy

- Designed to reduce colic and discomfort

### Finding the right nipple is important

- Finding the right nipple is important

# Highlights

## Anti-colic valve



Anti-colic valve is designed to keep air away from baby's tummy during feeding, to help reduce colic and discomfort.

## Natural latch on



The wide, soft and flexible nipple is designed to mimic the shape and feel of a breast, helping baby to latch on and feed comfortably.

## Finding the right nipple is important



If your newborn baby is consistently not taking enough milk throughout the feeding sessions or has complications in getting milk, switch to a nipple with a higher flow rate. If persistent feeding issues occur, consult a healthcare professional.

## Natural Response Nipple



The Natural Response Nipple works with your baby's natural feeding rhythm, making it easy to combine breast feeding and bottle feeding. The nipple has a unique opening which only releases milk when the baby actively drinks. So when they pause to swallow and breathe, the milk pauses too.

## Find the right nipple flow



Every baby feeds differently, and develops at their own pace. We've designed a range of flow rates so you can find the perfect one for your baby and personalize your bottle. All of the Natural Response Nipples are made of soft silicone.

## No-drip nipple



The nipple opening is designed to release milk only when the baby is feeding. So you can confidently avoid milk loss whether at home or on the go.

## Compatible across the range




Mix and match our breast pump, bottle and cup parts, and create the product that works for you, when you need it.

## Easy to hold



The ergonomic bottle is easy to grip at any angle to give maximum comfort during feeding. Easy for your hands and tiny hands to hold.

## Simple to use and clean



Wide bottle neck makes filling and cleaning easy. Only a few parts for quick and simple assembly.

# Specifications

## Material

**Bottle:** Polypropylene, BPA free\*

**Nipple:** Silicone, BPA free\*

## What is included

**9oz/260ml Baby Bottle:** 4 pcs

## Functions

**Nipple Features:** Natural latch on, No-drip design, Soft and flexible, Anti-colic valve

**Bottle design:** Wide neck, Ergonomic shape

**Bottle ease of use:** Easy to assemble, Easy to clean, Dishwasher safe, Easy to hold

## Development stages

**Stage:** 0-12 months

