

# PHILIPS

## sonicare



Oral Irrigator

Cordless Power  
Flosser 2000



HX3822/11

## Removes up to 99.9% of plaque in treated areas\*

### A thorough clean in-between

Effortlessly thorough cleaning between teeth to improve your gum health. Adjustable intensity and 360° rotation adapts to your needs to clean away plaque and food debris in hard-to-reach areas.

#### Easy and convenient

- Easy-to-fill 250 ml reservoir
- Fast universal charging with USB-A to Small Plug cable
- Cordless, compact design makes flossing simple
- Powerful and gentle for a thorough clean

#### Adapts to your needs

- 3 intensity settings to suit your needs
- 360° rotation for hard-to-reach spots


# Cordless Power Flosser 2000

Oral Irrigator

HX3822/11


## Highlights

### A powerful, gentle clean




Designed to remove plaque and food debris even from hard-to-reach areas. Leave your mouth feeling fresh and clean without the hassle of string flossing.

### Cordless, compact design



Compact, ergonomic design offers a convenient solution for everyday flossing. Cordless handle allows for easy portability and flossing anywhere.

### 360° nozzle rotation




Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation.

### 3 intensity settings



Choose the clean that suits your smile. Adjustable intensity settings adapt for comfort and cleaning power.

### 250ml reservoir



Reservoir capacity of 250 ml holds enough water for a recommended 60-second clean. Simply twist and remove or use the side-fill door for easy refilling.

### Fast charging with USB-A to Small Plug cable



Convenient charging with USB-A to Small Plug cable. You can clean for up to 14 days on just one charge.

## Specifications

### Consumer Trade Item

Height: 24.7 cm  
Width: 13.4 cm  
Length: 9.2 cm  
Net weight: 0.386 kg  
Gross weight: 0.604 kg  
GTIN: 8720689011051  
Country of origin: CN  
Harmonized system code: 96032100

### Outer Carton

Height: 27.2 cm  
Width: 15.1 cm  
Length: 29.9 cm  
Net weight: 1.158 kg  
Gross weight: 2.112 kg  
GTIN: 18720689011058

### Pallet

Height: 96.6 cm  
Width: 80 cm  
Length: 120 cm  
Net weight: 69.48 kg  
Gross weight: 151.98 kg  
GTIN: 68720689011053

### Service

Warranty: 2-year limited warranty

### Ease of use

Handle: Slim ergonomic design  
Nozzle attachment: Easily clicks on and off

### Performance

Cleaning: Clean entire mouth in 60-90sec

### Items included

Cordless Power Flosser: 1  
F1 Standard nozzle: 1  
USB charging cable: 1

### Modes

Clean: For a thorough everyday clean  
Intensities: 3

### Health benefits

Plaque removal: Removes up to 99.9% of plaque in treated areas\*

### Design and finishing

Color: White

