

PHILIPS

IPL hair removal
system


Lumea Precision

For use on body and face

15 minutes to treat lower legs

Lifetime >140.000 light
pulses

Cordless design

 Product Picture

SC2006/11



The ideal way to prevent hair regrowth

for smooth skin at home

Philips Lumea IPL is the most effective cordless IPL solution and works wonders to prevent the reappearance of hair on body. Gentle pulses of light, applied regularly, keep skin silky-smooth every day.

Long lasting smoothness

- Enjoy smooth skin every day
- Effortlessly effective

Hair removal with IPL technology

- Adapted for safe use in the comfort of your own home
- Suitable for a wide variety of hair and skin types

Gentle and effective application for face and body

- Precision attachment for safe facial treatment
- Bigger body attachment for fast treatment
- Slide and Flash mode for a more convenient application
- Safe and effective even on sensitive areas
- Five adjustable light energy settings

Maximum freedom, minimum hassle

- Cordless convenience
- No replacement parts, no hidden costs

Highlights

Long lasting smoothness



Philips Lumea applies gentle pulses of light to the hair root. As a consequence the hair sheds naturally and hair regrowth is inhibited. Repeating this treatment regularly will keep your skin smooth every day.

Effortlessly effective



Our clinical studies have shown significant hair reduction in as little as four bi-weekly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth.

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system.

For more than 10 years we carried out extensive consumer research with more than 2000 women.

Suitable hair and skin types

Suitable hair and skin tones

Skin Tone	White	Beige	Light brown	Medium brown	Dark brown	Brownish black and darker
Hair Color						
Black	✓	✓	✓	✓	✓	✗
Dark brown	✓	✓	✓	✓	✓	✗
Brown	✓	✓	✓	✓	✓	✗
Dark blond	✓	✓	✓	✓	✓	✗
Light blond	✗	✗	✗	✗	✗	✗
Red	✗	✗	✗	✗	✗	✗
White/grey	✗	✗	✗	✗	✗	✗

● This product is suitable for you to use ☹ This product is neither suitable nor safe for your skin and hair type

Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blond or white/grey hair. Philips Lumea is also NOT suitable for very dark skin.

Five light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.

Application areas



Use with confidence to remove hair on face (upper lip, chin and sideburns) and body parts including legs, underarm, bikini area, belly and arms.

Cordless convenience



Lumea's unique cordless design is perfect for maximum reach and ease of use.

No replacement parts needed



No replacement lamps or gel needed

Philips Lumea is a complete solution that comes ready to work right out of the box. More importantly, it does not require any replacement lamps or gels.

Extra safety for the face



The precision attachment is equipped with an extra integrated light filter and enables a safe and convenient application on the upper lip, chin and sideburns. It can also be used on other sensitive areas on the body, such as the underarms, the bikini area and hard-to-reach areas.

Specifications

Application areas

Body areas: Legs, Arms, Belly, Bikini, Armpits

Face areas: Chin, Sideburns, Upper lip

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light

Safety ring system: Prevent unintentional flashing

5 light energy settings: Adjustable to your skin type

Application mode

Slide and Flash: For convenient application

Step and Flash: For treatment on curvy areas

Corded / cordless use: Cordless

Power

Voltage: 100-240 V

Charging time: 100 min

Running time: Minimum 370 flashes

Technical specifications

High performance lamp: Delivers <gt/> 140,000 flashes

Attachments

Precision attachment (2cm²): For safe application on face

Body attachment (4cm²): For use below the neckline

Service

Warranty: 2 years global warranty

Items included

Charger: Multi voltage charger

Instructions for use: Quick Start Guide, User manual

Storage: Pouch

Application time

Lower leg: 8 min

Armpit: 1 min

Bikini line: 1 min

Face areas: 1 min

