

PHILIPS

sonicare

Standard sonic
toothbrush heads

S Sensitive

3-pack

Standard size

Click-on

For sensitive teeth and gums

HX6053/63



Ultra-sensitive cleaning. Ultra-soft comfort.

Superior cleaning for sensitive teeth and gums

The Philips Sonicare S Sensitive toothbrush head is perfect for those who seek relief from painfully sensitive teeth and gums. You'll get an effective clean that's easy on your gums.

An easy addition to your oral care routine

- Get the most effective plaque removal, every day

Innovative technology

- Philips Sonicare's advanced sonic technology
- Innovative bristle design softly sweeps plaque away

Designed to optimize performance

- Click-on design for simple brush head placement

An unrivaled deep clean

- Removes more plaque than a manual toothbrush


Highlights

Say goodbye to plaque



Densely-packed, high-quality bristles are clinically proven to remove more plaque than a manual toothbrush

Softly sweeps away plaque



This Philips Sonicare sensitive brush head features ultra-soft bristles for gentle, effective cleaning. And the contoured bristles fit the shape of your teeth, for an ultra-sensitive brushing experience. Also available in smaller, compact size for precision cleaning.

Easy to click-on



Your Sensitive brush head fits perfectly with any Philips Sonicare toothbrush handle, except PowerUp Battery and Essence. Simply click on and off for easy replacement and cleaning.

Reminder bristles



Brush heads become less effective after 3 months of use, but we'll remind you before this happens. Just monitor your blue replacement bristles, and when they turn white you'll know it's time for a fresh brush head.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Specifications

Design and finishing

Bristle stiffness feel: Ultra soft

Color: White

Reminder bristles: Blue bristle color fade away

Size: Standard

Compatibility

Brush head system: Click-on

Suitable for these models: 2 Series plaque control, 3 Series gum health, DiamondClean, EasyClean, FlexCare, FlexCare Platinum, FlexCare+, for Kids, HealthyWhite, HealthyWhite+

Items included

Brush heads: 3 S Sensitive standard

Quality and performance

Replacement: Every 3 months

Tested: for optimal usage

Health benefits

Plaque removal: Helps remove plaque

