

# PHILIPS

## Neck Massager

EMS microcurrent massage

2 hot compress levels

Flexible support structure

PPM3311



## Double care for the neck

### Support to soothe and massage to relax

Designed to relieve neck pressure, it combines micro current massage and warm compress. Anti-bow pattern can strongly support the chin, with neck support. Various modes and intensities are available, which is suitable for different users.

#### Combine head & neck support and massage

- Breakthrough combination of head & neck support and massage
- Anti-bow V-shaped design
- Flexible support structure
- Soft, highly elastic foam

#### Multiple massage method

- EMS dual frequency pulse massage
- 2 hot compress levels
- 4 major modes and 12 levels of intensity

#### Worry-free use

- Soft neck extension strap
- 150g lightweight design
- 120 minutes of long-lasting battery life

# Specifications

## General specifications

Massage type: EMS, Hot compress

Product Dimension (LxWxH): 170x180x110 mm

Net weight: about 0.15kg

Colour: Dark blue

## Main parameters

Auto shut-off time: 10 minutes

Working voltage: 3.7 V

Working power: 3 W

Battery: 950 mAh

Massage modes: 4

Charging current: 1 A

Charging port: USB Type-C

Charging voltage: 5 V

Massage strength levels: 12

## Service

Warranty period: 2 years

## Country of origin

Made in: China

