

Philips  
Neck Massager

EMS microcurrent massage  
2 hot compress levels  
Flexible support structure

PPM3311



## Double care for the neck

Support to soothe and massage to relax

Designed to relieve neck pressure, it combines micro current massage and warm compress. Anti-bow pattern can strongly support the chin, with neck support. Various modes and intensities are available, which is suitable for different users.

### Combine head & neck support and massage

- Breakthrough combination of head & neck support and massage
- Anti-bow V-shaped design
- Flexible support structure
- Soft, highly elastic foam

### Multiple massage method

- EMS dual frequency pulse massage
- 2 hot compress levels
- 4 major modes and 12 levels of intensity

### Worry-free use

- Soft neck extension strap
- 150g lightweight design
- 120 minutes of long-lasting battery life

**PHILIPS**

## Neck Massager

EMS microcurrent massage 2 hot compress levels, Flexible support structure

# Specifications

### General specifications

- Product Dimension (LxWxH): 170x180x110 mm
- Net weight: about 0.15kg
- Massage type: EMS, Hot compress
- Colour: Dark blue

### Main parameters

- Working voltage: 3.7 V
- Auto shut-off time: 10 minutes
- Working power: 3 W
- Massage modes: 4

- Battery: 950 mAh
- Charging current: 1 A
- Charging port: USB Type-C
- Charging voltage: 5 V
- Massage strength levels: 12

### Country of origin

- Made in: China

### Service

- Warranty period: 2 years

PPM3311DB/90

# Highlights

### Combination of head and neck support and massage

The newly created V-shaped anti-duck structure incorporates EMS massage and hot compress. While massaging, it helps you straighten your cervical spine and sit upright, creating a multi-functional massage experience.

### Anti-bow pattern

The area close to the chin adopts a newly created V-shape, which is ergonomically designed and fits the curve of the chin. Support your chin steadily to prevent excessive bowing and reduce pressure on your head and neck.

### Flexible support structure

The built-in curved elastic support rack can give the head a good assistance after wearing this motorcycle. Hold the chin to prevent excessive bending and relieve stress.

### Soft, highly elastic foam

The inner core of slow rebound foam has just the right elasticity and is suitable for different neck lengths. The foam has multiple breathable holes, and the soft and smooth fabric on the device makes your massage relaxing, refreshing and not stuffy.

### EMS microcurrent massage

The medium and low dual-frequency pulses are output uniformly and stably, simulating various massage techniques such as pressing, tapping, kneading, and massage, and mobilizing neck muscle contraction to relax.

### 2 hot compress levels

Embedded with a thin heating sheet, gentle heat penetrates into the core muscle area of the neck. You can enjoy hot compress while massaging to relax the tight neck.

### 4 major modes and 12 levels of intensity

Different modes and pulse intensity levels can be freely combined. The intensity can be strong or mild, giving you a diverse massage experience. The buttons are designed with concave and convex patterns for easy and convenient operation.

### Neck extension strap

The neck strap is lengthened and can be attached in two ways to comfortably fit a variety of neck sizes to meet the needs of different people.

### Lightweight design

The whole device weighs only 150g. It is light and easy to wear. It does not put additional burden on the cervical spine and brings a burden-free massage experience.

### Long lasting battery life

It has a built-in 950mAh lithium battery, which can last for 120 minutes on a full charge. It is equipped with a Type-C charging cable to recharge the battery at any time, making it a good travel companion.



Issue date 2024-01-19

Version: 2.2.2

EAN: 48 97126 05095 0

© 2024 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)