



Philips Axon



Philips Axon is the next generation of reliable, serial-to-network bridges within the Philips Capsule Medical Device Information Platform (MDIP), enabling the automatic collection of medical device data. The Philips Axon features an always-on design and connects to the hospital network through a wired Ethernet or standard 802.11 a/b/d/g/h/j/n/ac/r/w wireless network.



Features

Flexible configurations to meet your needs

- One, four, or eight port configurations to fit a wide range of use cases
- Support for both wired and wireless network connections to support stationary or mobile care use cases
- Small profile design enabling deployment in crowded environments

Easy and simple hardware deployment

- Power over ethernet offers the possibility for a single cable to support both power and network connection
- True medical device plug and play integration

Reliable connectivity and data integration

- Enables the integration of device data into the hospital EMR and other downstream systems
- Brings “always-on” data feeds to support patient surveillance capabilities to care areas with minimal space
- Expands access to data-driven interventions and precision care protocols in more areas of the hospital

For more information, contact us

North America

+1 800-260-9537

support@capsuletech.com

International Offices

+33 1 84 17 12 50

international@capsuletech.com

CapsuleTech.com

