

PHILIPS

sonicare

Power toothbrush

For Kids Design a
Pet Edition

Sonicare technology

Customisable stickers

Free Sonicare for Kids app



HX3603/01



Helps kids prevent cavities

A fun way to build lifetime, healthy habits

Help build healthy brushing habits with Sonicare for Kids. It's safe and gentle on growing teeth and gums and helps prevent cavities. Our special Design a Pet edition comes with fun stickers for kids to make it their own.

Designed around you

- Makes traveling easier for you

Make brushing fun

- Create your own personal design with fun, reusable stickers
- Free Sonicare for Kids app coaches them to brush better

Helps kids prevent cavities

- Sonicare technology helps prevent cavities

Build lifetime healthy habits

- Brush heads are specially designed for kids
- SmarTimer and KidPacer keep kids brushing for 2 minutes
- 2 intensity settings for a gentle clean

Highlights

Sonicare technology



Gentle pulsing action reaches deep between teeth and along the gum line. So, kids get the best possible clean while they are still learning to brush.

Brush head designed for kids



The brush heads are specially designed to help kids get a thorough clean while they learn to brush. One compact brush head is included. Standard size heads sold separately.

Create your own personal pet



Special Design a Pet edition includes reusable stickers so that kids can turn their toothbrush into new animal characters any time. Let's make brushing a fun part of their day!

SmarTimer and KidPacer




Coach kids to build healthy habits. SmarTimer keeps them brushing for a full 2 minutes and KidPacer makes sure that they switch to a new area of the mouth every 30 seconds.

2 intensity settings



Developing teeth need extra protection. Choose between Gentle and Extra Gentle settings, with a rubberised brush head for a comfortable clean.

Free Sonicare for Kids app



Build healthy habits to last a lifetime. Our free Sonicare for Kids app supports kids as they brush twice a day for the full dentist-recommended 2 minutes and have fun doing it!

Slim travel case



Our premium travel case lets you store your toothbrush hygienically, while our compact USB charging base keeps you topped up when you are on the go. You get two weeks of regular use from a single full charge, but the charger's just the ticket for longer trips.

Specifications

Power

Voltage: DC5V

Technical specifications

Operating time (full to empty): 14 days

Battery: Rechargeable

Battery type: Lithium ION

Consumer Trade Item

Height: 22

Width: 13.9

Length: 6.3

Net Weight: 0.291

Gross Weight: 0.4

GTIN: 8720689018777

Country of origin: CN

Harmonized system code: 85098000

Outer Carton

Height: 25.4

Width: 20.8

Length: 15.5

Net Weight: 0.873

Gross Weight: 1.363

GTIN: 18720689018774

Design and finishing

Colour: White with Aqua Blue button

Service

Warranty: 2-year limited warranty

Other

Height: 115

Width: 80

Length: 120

GTIN: 68720689018779

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

Items included

Brush heads: 1 Sonicare for Kids compact

Travel case: 1 slim travel case

Charger: 1 USB charger

Handle: 1 Sonicare for Kids Pet handle

Stickers: 2 pet themed sticker sheets

Cleaning performance

Timer: SmarTimer and KidPacer

