



Vue PACS 12.2.8

IHE Integration Statement

Part # HA1668
2021-03-26

Integrating the Healthcare Enterprise (IHE) Integration Statement			
Vendor	Product Name	Version	Date
Philips Medical Systems	Vue Archive (Vue for Clinical Collaboration Platform), Vue PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry.	Vue PACS, Vue Archive v12.2.8	2021-03-10
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed in this table.			

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None
Imaging Object Change Management (IOCM)	Image Manager/Image Archive Change Requestor	None
Access to Radiology Information (ARI)	Image Manager/Image Archive Image Display Report Repository External Report Repository Access	None
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None
Cross-enterprise Document Sharing (XDS.b)	Document Consumer Document Repository Document Registry Document Source	Multiple document
Cross-enterprise Document Sharing for Imaging (XDS.b-1)	Imaging Document Consumer Imaging Document Source	Set of DICOM Instances, Multiple documents None
Cross-enterprise User Assertion (XUA)	X-Service User	None
Cross-enterprise Document Reliable Interchange (XDR)	Document Recipient Document Source	Accept Limited Metadata
Cross-Enterprise Sharing of Scanned Documents (XDS-SD)	Content Consumer	None
Basic Patient Privacy Consents (BPPC)	Content Consumer	None

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Demographic Query HL7 V3 (PDQv3)	Patient Demographics Consumer	None
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None
Evidence Documents (ED)	Image Display Image Manager/Image Archive	None
Cardiology Evidence Documents (ED-CARD)	Image Display Image Manager/Image Archive	None
Simple Image and Numeric Report (SINR)	Report Repository External Report Repository Access	None
Mammography Image (MAMMO)	Image Manager/Image Archive Image Display	None
Key Image Note (KIN)	Evidence Creator Image Display Image Manager/Image Archive	None
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None
Audit Trail and Node Authentication (ATNA)	Secured Application Secured Node	None
Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web Content None None None
Reporting Workflow Profile (RWF)	Image Manager/Image Archive	None
Import Reconciliation Workflow Integration Profile (IRWF)	Image Manager/Image Archive	None
Teaching File and Clinical Trial Export (TCE)	Receiver	None
Post-Processing Workflow (PWF)	Image Manager/Image Archive	None
Presentation of Grouped Procedures Integration Profile (PGP)	Image Manager/Image Archive	None

IHE Conformance Statements	Philips IHE (www.philips.com/IHE)
DICOM Conformance Statements	Philips DICOM (www.philips.com/DICOM)
More about products from Philips Healthcare	https://www.philips.com/healthcare/about/customer-support
General information on IHE	<ul style="list-style-type: none"> • http://www.ihe.net/ (International and North America) • https://www.ihe-europe.net/ (Europe) • http://www.ihe-j.org/en/ (Japan)



Philips Medical Systems Nederland B.V.
Veenpluis 6, 5684 PC Best, The Netherlands



© 2021 Koninklijke Philips N.V.

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Philips or its licensors. No title or ownership in the Software is conferred to you. Use of the Software is subject to the end user license conditions as are available on request.

HA1668_A/ * 2021-03-26 – en-US

End of Document