

# PHILIPS

## AVENT

### Soothie Shapes pacifier

One piece silicone design

0-3m

Orthodontic & BPA-Free

2-pack



SCF194/02



## Promotes natural suckling and bonding

### Flexible medical grade silicon

The Philips Avent Soothie Shapes pacifier gives you an additional way to bond with your little one as our unique one piece design lets you give the comfort of your finger. Made of medical grade silicon, Soothie Shapes is easy to clean.

#### Easy to clean

- Easy to sterilize for extra hygiene

#### Comfortable round nipple

- Designed for natural oral development

#### Fingerheld solution

- Unique design supports bonding

#### Medical grade silicone

- Made of BPA-free, durable, medical grade silicone

#### Used by medical professionals

- Distributed in hospitals in the USA


# Highlights

## Used in hospitals



Trusted by healthcare professionals to calm newborns, Soothies are distributed in hospitals in the USA.\*

## Unique shape



Soothie Shapes is a little different to our other pacifiers. Its unique shape lets you place your finger in the nipple so you can bond with your baby by helping them suckle.

## Orthodontic nipple



Our collapsible silicone nipple has a symmetrical shape that respects your baby's palate, teeth and gums as he or she grows.

## Medical grade silicone



The pacifier is made of BPA-free, durable, yet flexible medical grade silicone.

## Easy to sterilize



Keeping your little one's Soothies and pacifiers clean is easy. Simply put them in your sterilizer or submerge in boiling water.

# Specifications

## Hygiene

Can be sterilized  
Easy to clean

## Consumer Trade Item

Height: 11.5 cm  
Width: 10.3 cm  
Length: 4.8 cm  
Net weight: .038 kg  
Gross weight: .068 kg  
GTIN: 8710103742999  
Country of origin: US  
Harmonized system code: 392690

## Outer Carton

Height: 12.9 cm  
Width: 11.6 cm  
Length: 32.6 cm  
Gross weight: .896 kg  
GTIN: 18710103742996

## Safety

BPA-free

## What is included

Soothie Shapes pacifier: 2 pcs

