

PHILIPS

Neck Massager

EMS microcurrent massage

2 hot compress levels

Flexible support structure

PPM3311



Double care for the neck

Support to soothe and massage to relax

Designed to relieve neck pressure, it combines micro current massage and warm compress. Anti-bow pattern can strongly support the chin, with neck support. Various modes and intensities are available, which is suitable for different users.

Combine head & neck support and massage

- Breakthrough combination of head & neck support and massage
- Anti-bow V-shaped design
- Flexible support structure
- Soft, highly elastic foam

Multiple massage method

- EMS dual frequency pulse massage
- 2 hot compress levels
- 4 major modes and 12 levels of intensity

Worry-free use

- Soft neck extension strap
- 150g lightweight design
- 120 minutes of long-lasting battery life

Specifications

General specifications

Massage type: EMS, Hot compress

Net weight: about 0.15kg

Product Dimension (LxWxH): 170x180x110 mm

Colour: Dark blue

Main parameters

Auto shut-off time: 10 minutes

Working voltage: 3.7 V

Working power: 3 W

Massage modes: 4

Battery: 950 mAh

Charging current: 1 A

Charging port: USB Type-C

Charging voltage: 5 V

Massage strength levels: 12

Service

Warranty period: 2 years

Country of origin

Made in: China

