




Philips
Infrared lamp

150 W

For targeted treatment
Adjustable angle



PR3110/00

Effective*, drug-free pain relief

Clinically proven technology*

Our clinically proven technology provides effective relief from muscle and joint pains without the use of drugs. It stimulates local blood circulation and increases the local supply of oxygen to provide relaxation to the injured area

Proven health benefits

- Clinically proven
- Classified medical device

Targeted treatment

- Localized pain relief

Effective* pain relief

- Pain relief with better mobility
- Helps the body to heal in a natural manner

PHILIPS

Infrared lamp
150 W For targeted treatment, Adjustable angle

Specifications

Zone treatment

- Treatment area: diameter of 20 cm
- PAR38 incandescent lamp: 150 W

Technical specifications

- Voltage: 220-240 V
- Cord length: 250 m
- Power: 150 W
- Frequency: 50-60 Hz
- Insulation: Class II (double isolation)
- Lifetime of lamp: 300 hour(s)

Weight and dimensions

- Product weight: 1 kg
- Product dimensions: 21,2 x 16,0 x 21,5 cm
- F-box dimensions: 16,3 x 22,4 x 22,0 cm

Logistic data

- Country of origin: Hungary

Safety

- IEC certified: IEC 60601-1, 3rd edition
- 30 minutes auto shut-off

Easy to use

- Adjustability: 0-40 degrees (backwards)
- 1 handgrip
- Mechanical timer

Medical appliance

- Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

PR3110/00

Highlights

Effective* pain relief

The InfraCare PR3110 provides an effective focused and targeted treatment by increasing the local blood circulation to mobilize your body against pain. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Focused infrared warmth

InfraCare enables the body to heal in a natural manner. The PR3110 InfraCare device helps to provide relief from pain by focused treatment, which helps in proven improved mobility.

Clinically proven

Trusted technology used by professionals. Clinically proven to objectively lead to functional joint improvement*. InfraCare has proven to lead to significant pain reduction in 71% of rheumatoid arthritis patients and in 72% of low back pain patients*.

Classified medical device

A safe and effective way to treat your pain. InfraCare complies with IEC 60601-1 3rd edition for medical devices.

Adjustable angle

The InfraCare 150W Lamp can be easily positioned from 0-40 degrees to relieve pain in a comfortable and an effective manner in small areas such as hips, arms and knees.



Issue date 2024-04-20

Version: 2.0.2

EAN: 87 10103 79744 9

© 2024 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* W. Siems, et al (2010): "Pain and mobility improvement and MDA plasma levels in degenerative osteoarthritis, low back pain, and rheumatoid arthritis after infrared A-irradiation", Acta Biochimica Polonica 57(3), 313-319