


PHILIPS

Imaging


ChromaFlo

Image interpretation made easy


ChromaFlo aids image interpretation



Decide



Guide



Confirm

Effortless to use

- Simply press the button
- No special procedures, no flushing during imaging
- Works with Philips' fast, plug-and-play digital IVUS catheters as shown below

Broad versatility

- Appropriate for peripheral vessels including bifurcations, SFA, iliac arteries and below-the-knee
- Provides both lumen and vessel information
 - Colorization highlights the lumen
 - Automatically co-registered grayscale IVUS allows observation of plaque burden, calcium, and vessel size

Compatible Philips digital IVUS catheters



Visions PV .014P RX
.014" guide wire compatible

Visions PV .018 RX
.018" guide wire compatible

Reconnaissance PV .018 OTW
.014" & .018" guide wire compatible

Always read the label and follow the directions for use.
Products subject to country availability. Please contact your local sales representative.

© 2024 Koninklijke Philips N.V. All rights reserved.
Trademarks are property of Koninklijke Philips N.V. or their respective owners.

Approved for external use. D2024052905 (D041343-02) * APR 2024

www.philips.com/IGTdevices



Volcano Corporation
3721 Valley Centre Drive, suite 500
San Diego, California, 92130, USA

Volcano Europe SA/NV
Excelsiorlaan 41
B-1930 Zaventem, Belgium

