

PHILIPS

sonicare



Dental Professional Trial Unit

7100

Up to 20 x more plaque removal

4 modes, 3 intensities

Visible pressure alert

Philips Sonicare app



HX7428/02

Perfect clean with up to 7 x healthier gums**

Up to 20 x more plaque removal*

The Philips Sonicare 7100 features Next-Generation Sonicare technology, 4 brushing modes, 3 intensity levels, a visible pressure alert and our G3 Premium Gum Care brush head.

Premium, gentle cleaning

- One brush head for complete care

Brushing made effortless

- Personalise your clean with 12 brushing settings
- Protect gums with the visible pressure sensor
- Always know when it's time for a new brush head
- Elevate patients' routine with the Philips Sonicare app

Superior clean

- Next-Generation Sonicare Technology

Designed for your lifestyle

- Laser-welded for strength and durability
- Make travel easy with our compact, ultra-durable travel case
- Up to 3 weeks of brushing with one full charge
- Optimise brushing with Smartimer and BrushPacer


Highlights

A3 All-in-One Brush Head +more



The A3 Premium All-in-One brush head is designed with patients in mind. With the A3 brush head, patients can achieve up to 20 x more plaque removal*, up to 15 x healthier gums** and remove up to 100% more stains***. Also included: The S2 Sensitive Brush Head which is our gentlest and softest brush head, and is clinically proven to provide a gentle clean that removes 10 x more plaque in hard-to-reach areas* The G3 Premium Gum Care Brush Head provides up to 7 x healthier gums in just 2 weeks****.

Next-Gen Sonicare technology




Enjoy a delightful brushing experience and new level of oral care with Next-Generation Sonicare technology. While you brush, the toothbrush's innovative magnetic system adjusts the brushing power during each cleaning session, delivering optimal performance of 62,000 bristle movements per minute, even in hard-to-reach areas. Meanwhile, dynamic fluid action supports the bristles to clean plaque, helping you achieve a complete clean every time you brush.

12 brushing settings



Patients' preferred brushing is in their hands with the Philips Sonicare 7100. Four brushing modes help take care of gum health, whitening, sensitivity and everyday cleaning. Gum Health mode supports healthy gums by gently massaging along the gumline. White mode works to remove surface stains and polish teeth for a whiter smile. Sensitive mode provides gentle, thorough cleaning, while Clean mode focuses on plaque removal for a superior clean. For a personalised clean, patients can combine any mode with one of three intensity levels, from High to Low, for a brushing session that's just right.

Light-ring pressure sensor




Patients may not notice when they're brushing too hard, but this toothbrush will. The smart optic sensor precisely detects excessive brushing pressure, then flashes purple while reducing vibration. Patients will always know when to adjust their brushing, and their gums will stay protected.

Philips Sonicare app



The Philips Sonicare app seamlessly pairs with the toothbrush and supports patients in improving and maintaining their oral health. Our app keeps track of brushing progress, providing personalised guidance and progress reports so patients can gain valuable insights about their brushing.

2-minute timer and brushing pacer



Two minutes is all patients need for a complete clean, so our Smarttimer lets patients know when time's up. For an even clean throughout your mouth, the BrushPacer divides brushing sessions into segments, with a slight pause in vibration being a sign to move on. This supports patients in achieving an even clean every time they brush.


BrushSync replacement reminder



Replacing your brush head regularly is an important part of maintaining good oral health. With BrushSync technology, it's easy for patients to keep their brush head performing at its best. By tracking usage frequency and brushing intensity, the innovative technology then reminds patients when it's time for a new brush head.

Highlights

Designed to last



We've re-designed our Philips Sonicare toothbrushes from the inside out to be strong, lightweight and long-lasting. Thanks to our advanced laser-welding process, our toothbrushes are designed to brush patients' teeth for years to come.

Up to 3 weeks between charges



This toothbrush uses a high-performance battery that is engineered to enhance the patient's brushing experience. With a single charge, patients can enjoy up to three weeks of regular brushing before needing to recharge it.

Compact travel case



Our redesigned travel case is compact enough to slip into a bag and strong enough to keep your Philips Sonicare toothbrush protected along the way. Not only does it keep the toothbrush protected, but it also conveniently holds two brush heads. Now patients can travel confidently wherever their journey takes them.

Specifications

Intensities

High: To boost your clean

Medium: For an everyday clean

Low: For sensitive teeth and gums

Product Update

Support Period: Philips will provide the necessary security updates until 31 December 2026.

Power

Voltage: 100-240 V

Technical specifications

Operating time (full to empty): 21 days

Battery: Rechargeable

Battery type: Lithium ION

Energy consumption: Standby without display
<0.5 W

Design and finishing

Colour: White

Service

Warranty: 2-year limited warranty

Compatibility

Bluetooth wireless technology: Connected brushing app

Android compatibility: Android phones with Bluetooth 4.0 and above

iOS compatibility: iPhone with Bluetooth 4.0 and above

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Battery indicator: Illuminated icon indicates battery life

Timer: SmarTimer and BrushPacer for complete and thorough brushing

USB-A Charging: Power adapter not included

Items included

Charger: 1

Travel case: 1 slim travel case

Handle: 1 7100 rechargeable toothbrush

Brush head: 1 A3 Premium All-in-One, 1 S2

Sensitive, 1 G3 Premium Gum Care

USB-C Cable: 1

Cleaning performance

Plaque removal: Up to 20 x more effective*

Speed: 62,000 brush movements/min

Gum Health: Up to 15 x healthier gums***

Whitening: Up to 100% fewer stains**

Modes

Clean: For exceptional everyday cleaning

White: To remove surface stains

Sensitive: For sensitive teeth and gums

Gum Health: Gently massages your gums

Smart sensor technology

Pressure feedback: Light ring illuminates purple, Reduced vibration and pulsing sound

Replacement Reminder: BrushSync lets you know when to replace your brush head



* vs. a manual toothbrush.

** in Clean mode vs. a manual toothbrush.

*** vs. a manual toothbrush. In-vitro study; actual results may vary.

**** in Gum Health mode vs a manual toothbrush; measured by GBI.