

PHILIPS

sonicare



Rechargeable Sonic Toothbrush

6100

Up to 20 x more plaque removal

2 modes, 3 intensities

Visible pressure alert

21-day battery life



HX7408/01

Gentle clean, even for sensitive teeth and gums

Up to 20 x more plaque removal*

The Philips Sonicare 6100 features Next-Generation Sonicare technology, a dedicated Sensitive mode, 3 brush heads for plaque, sensitivity and whitening, and a low intensity. Plus, a pressure sensor for when patients brush too hard.

Gentle, superior cleaning

- Complete care for a healthier mouth

Brushing made effortless

- Protect gums with the visible pressure sensor
- Always know when it's time for a new brush head
- Personalise your clean with 6 brushing settings

Superior clean

- Next-Generation Sonicare Technology

Designed for your lifestyle

- Up to 3 weeks of brushing with one full charge
- Laser-welded for strength and durability
- Optimise brushing with Smartimer and BrushPacer

Highlights

A3 All-in-One and S2 Sensitive



The A3 Premium All-in-One brush head is designed with patients in mind. With the A3 brush head, patients can achieve up to 20 x more plaque removal*, up to 15 x healthier gums** and remove up to 100% more stains***. Also included: Thanks to the new S2 Sensitive Brush Head made using Sonicare's softest bristles, patients can achieve a gentle clean that removes up to 10 x more plaque in hard-to-reach areas*.

Next-Gen Sonicare technology



Enjoy a delightful brushing experience and new level of oral care with Next-Generation Sonicare technology. While you brush, the toothbrush's innovative magnetic system adjusts the brushing power during each cleaning session, delivering optimal performance of 62,000 bristle movements per minute, even in hard-to-reach areas. Meanwhile, dynamic fluid action supports the bristles to clean plaque, helping you achieve a complete clean every time you brush.


6 brushing settings



Designed for sensitive teeth and gums, the Philips Sonicare 6100 allows patients to choose between two brushing modes. Sensitive mode

provides gentle, thorough cleaning, while Clean mode focuses on plaque removal for a superior clean. For a personalised clean, patients can combine any mode with one of three intensity levels, from High to Low, for a brushing session that's just right.

Light-ring pressure sensor




Patients may not notice when they're brushing too hard, but this toothbrush will. The smart optic sensor precisely detects excessive brushing pressure, then flashes purple while reducing vibration. Patients will always know when to adjust their brushing, and their gums will stay protected.

BrushSync replacement reminder




Replacing your brush head regularly is an important part of maintaining good oral health. With BrushSync technology, it's easy for patients to keep their brush head performing at its best. By tracking usage frequency and brushing intensity, the innovative technology then reminds patients when it's time for a new brush head.

2-minute timer and brushing pacer



Two minutes is all patients need for a complete clean, so our Smartimer lets patients know when time's up. For an even clean throughout the mouth, the BrushPacer divides brushing sessions into segments, with a slight pause in vibration being a sign to move on. This supports patients in achieving an even clean every time they brush.

Designed to last



We've re-designed our Philips Sonicare toothbrushes from the inside out to be strong, lightweight and long-lasting. Thanks to our advanced laser-welding process, our toothbrushes are designed to brush patients' teeth for years to come.

Up to 3 weeks between charges



This toothbrush uses a high-performance battery that is engineered to enhance your brushing experience. With a single charge, patients can enjoy up to three weeks of regular brushing before needing to recharge it.

Specifications

Intensities

High: To boost your clean

Medium: For an everyday clean

Low: For sensitive teeth and gums

Power

Voltage: 100-240 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): 21 days

Energy consumption: Standby without display
<0.5 W

Battery type: Lithium ION

Design and finishing

Colour: White

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Battery indicator: Illuminated icon indicates battery life

Handle compatibility: Easy click-on brush heads
Timer: SmarTimer and BrushPacer for complete and thorough brushing

USB-A Charging: Power adapter not included

Items included

Handle: 1 6100 rechargeable toothbrush

Brush head: 30 S2 Sensitive, 30 A3 Premium All-in-One

Charger: 1

Cleaning performance

Plaque removal: Up to 20 x more effective*

Gum Health: Up to 15 x healthier gums***

Whitening: Up to 100% fewer stains**

Speed: 62,000 brush movements/min

Modes

Clean: For exceptional everyday cleaning

Sensitive: For sensitive teeth and gums

Smart sensor technology

Pressure feedback: Light ring illuminates purple, Reduced vibration and pulsing sound

Replacement Reminder: BrushSync lets you know when to replace your brush head

