

PHILIPS

Mask

InnoSpire Go

Adult large

Soft seal mask

For use with InnoSpire Go

Replacement part

HH1344

 Product Picture

Aids in treatment

The detachable LiteTouch facemask has a soft contoured cushion that fits easily and rests gently on the face to provide a more comfortable treatment.*

Specifications

Package contents

Include: Large LiteTouch Mask

Material

BPA: Not manufactured with BPA

Product details

Packaging size: 82mm x 97mm x 100mm

Packaged weight: 80g

LiteTouch mask: For use with InnoSpire Go

Materials: Polyethylene terephthalate (PET),

Silicone

Size: Adult large

Warranty: 1 year

