

PHILIPS

sonicare



Sonic electric
toothbrush with app

DiamondClean
Smart 9300

Up to 20x more plaque removal*

Up to 15x healthier gums**

Up to 100% less stains***

Visible pressure alert



HX9903/25

Complete care for a healthier mouth

Feel confident with 100% coverage

The DiamondClean Smart provides a deep clean for a whiter, healthier smile. High-performance brush heads let you focus on all areas of your oral health and our Smart Sensor technology gives you personalized feedback and coaching.

Personalized in-app coaching to perfect brushing

- Smart guidance for gentle care
- Helps you improve your brushing technique
- Guides you and provides instant feedback as you brush
- Takes you back to any missed spots for 100% coverage
- Start taking better care of your trouble areas

Optimize your brushing

- Make sure you're getting the most from your brush head
- Automatically selects the best mode for you
- Premium charging glass and travel case
- Choose from 4 modes and 3 intensity settings

Premium performance

- All-in-One brush head for complete care
- Up to 4x more surface contact**** for an effortless clean
- Sonic brushing: the masterclass in effective cleaning

DiamondClean Smart 9300

Sonic electric toothbrush with app

HX9903/25

Specifications

Connectivity

Bluetooth® wireless technology: Connected brushing app

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): 14 days*****

Battery: Rechargeable

Battery type: Lithium ION

Design and finishing

Color: Pink

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones with Bluetooth 4.0 and above, Search Sonicare on Google Play Store for detailed OS compatibility.

iOS compatibility: iPhone with Bluetooth 4.0 and above, Search Sonicare on App Store for detailed iOS compatibility.

Ease of use

Display: Illuminated display

Battery indicator: Illuminated icon indicates battery life

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Brush head recognition: Syncs with optimal mode

Replacement reminder: To always ensure best results, reminder icon lights up

Items included

Handle: 1 DiamondClean Smart

Travel case: Premium travel case

Brush heads: 1 A3 Premium All-in-One, 1 C3

Premium Plaque Control

Glass charger: 1

Cleaning performance

Speed: Up to 62000 brush movement/min

Health benefits: Up to 15x healthier gums**

Performance: Removes up to 20x more plaque*

Timer: BrushPacer and SmarTimer

Whitening benefits: Up to 100% less stains***

Pressure feedback: Ring lights up in purple, Vibrates handle to alert user

TouchUp: Ensures 100% coverage

Modes

3 intensities: High, Medium, Low

Clean: For exceptional everyday clean

Deep Clean+: For an invigorating deep clean

Gum Health: Special attention to molars

White+: To remove surface stains

BrushSync Mode Pairing

C3 Premium Plaque Control: Pairs with Clean mode

A3 Premium All-in-One: Pairs with Clean mode

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

Location sensor: Tracks & improves coverage

Scrubbing sensor: Guides to reduce scrubbing


3D mouth map: Highlights trouble areas

TouchUp: To address missed spots

Position sensor: Guides to areas you missed


Highlights

TouchUp feature



If you happen to miss a spot while brushing your teeth, your app's TouchUp feature will show you. You can then go back for a second pass, and be confident you are getting a complete clean, every time.

Focus areas and goal-setting



Any trouble areas your dentist has pointed out? We'll highlight them on your in-app 3D mouth map, and remind you to pay extra attention to these areas.

Design with storage in mind



Charging a toothbrush has never been so effortless. Simply place your toothbrush into the sleek glass holder after brushing, and it instantly starts charging. Not only does your glass charge fit stylishly into your bathroom, it also doubles as a rinsing glass. Enjoy two weeks of regular use from a single full charge. Our premium travel case lets you store your toothbrush hygienically when you're on the go.



* vs. a manual toothbrush.

** when used with A3 brush head vs. a manual toothbrush.

***vs. a manual toothbrush. In a lab study, actual results may vary.

****than a DiamondClean brush head.

*****based on two periods of two-minute brushings per day on standard mode