



Lumify**



Compact 5000 series



Affiniti



EPIQ

Talk	●	●	●	●
Text		●	●	●
Screen share	●	●	●	●
Remote control *†		●	●	●
Remote diagnostic quality*			●	●
Asset sharing†		●	●	●
Multi-party*			●	●
System-to-system connection*			●	●
Remote user operating system				
– Windows (via Application, Chrome and Edge browsers)	●	●	●	●
– Mac OS with Chrome and Edge browsers*	●		●	●
– iOS*	●		●	●
– Android*	●		●	●

* Requires release 7.0.5 or higher. Diagnostic use and remote access via mobile device or browser requires release 9.0 or higher.

Multi-party and system-to-system connect require release 10.0 or higher.

** Lumify used with Android device.

† Available with Windows application only.