

PHILIPS

AVENT

Manual Breast Pump

Essential Manual
Breast Pump

Ergonomic expression handle

Lightweight design

BPA free

SCF417/11



Stimulate milk flow, express easily

For gentle and comfortable pumping

The Philips Avent Essential Manual breast pump designed to stimulate milk flow to help you express as much as milk possible. It helps to build back-up supply of breast milk so you can share the care to simply enjoy a moment of independence.

Easy to clean

- Easy to clean and set-up (dust cap)

Gentle and comfortable

- Soft and adaptive silicone cushion for gentle stimulation

Comfortable expression

- Ergonomic handle design for easy expression

Portable and convenient

- Lightweight and compact design

BPA free

- BPA free*


Essential Manual Breast Pump

Manual Breast Pump

SCF417/11


Highlights

Soft and adaptive silicone cushion



Soft and comfortable petal-shaped silicone massage cushion; gently stimulates milk flow while massaging the breast. It brings together gentle compression and suction to trigger fast let down.

Lightweight and compact design



Pump easily anywhere, any time with its compact and lightweight design, you can carry it in your bag anywhere you go. You can directly pump into Avent feeding bottles as the expression kit is compatible.

Ergonomic handle design for easy expression



With easy expression handle designed ergonomically, you can pump at your own pace and decide on your own rhythm simply by adjusting the vacuum level by pressing the handle at your own intensity.

Easy to clean and set-up



Designed for effortless use, featuring easy assembly, disassembly, and cleaning with its minimal, user-friendly parts. The dust cap will help to keep your breast pump clean and free from dust and dirt.

BPA free*



Made of durable and safe materials with 0% BPA.*

Specifications

Material

BPA free*

What is included

Breast milk container (125 ml/ 4 oz): 1 pcs

Expression kit with cushion: 1 pcs

Expression handle: 1 pcs

Dust cap: 1 pcs

Ease of use

Breast pump use: Easy cleaning & assembly, with lightweight design

