

# PHILIPS

Ultrasound

Affiniti Elevate

# What's new in Affiniti Elevate



Affiniti Elevate

Affiniti Elevate offers a combination of robust performance and value that is designed for your everyday.

**Welcome to the future of ultrasound.**

- 1 Higher diagnostic confidence**  
Now available on Affiniti Elevate, see **53%\* improved penetration** in superficial applications\*\* with the **mL26-8 transducer**
- 2 Extension of Next Gen Auto Scan**  
**Auto Scan** is now featured across **44 tissue-specific presets (TSPs)** for Affiniti Elevate
- 3 Extension of Flow Viewer**  
**Flow Viewer** is now featured across **74 TSPs** for Affiniti Elevate
- 4 Improved workflow**
  - Body marker imports
  - Color-coded annotations
  - Fusion and Navigation personalized touchscreen
  - Quick launch presets enhance user workflow

Now available on Affiniti Elevate, the **mL26-8 transducer improves axial resolution 55%\*** in superficial applications\*\*

**Quick launch presets**  
**Offers 50% decrease in the image optimization time** during an abdominal exam†

\*D001832939

\*\*Compared to the predecessor transducer L15-7io.

†D001833994, Marketing Claim Evidence for VM12.0 Workflow Efficiency Quick Launch Preset.

# Step into the future

Now your ultrasound system does more for you so you can keep enhancing patient care.



## Stay ahead with AI in shared service cardiovascular scanning

Everyday echo with **increased efficiency and clinical confidence**

AI-enabled Auto Strain and Auto Measure offer fast, reproducible results for easy comparison of functional data



## Speed diagnosis and surgical planning

FlexVue with measurements allows a **2D image to be easily extracted** from a 3D volume

Projected view of curved structures displayed in a single image



## Get connected for fast access to capabilities

**Remote Software Delivery (RSDi)** allows you to receive diagnostics and software remotely, schedule updates on your own time without system downtime and receive sustaining updates to keep your capabilities and security current.

# Ultrasound that shows you more, more easily

Stunning images, acquired with ease. New levels of clinical information and reproducibility. Versatility, with enhanced workflow for fast exams with diagnostic confidence. That's the future of ultrasound with Affiniti Elevate.



© 2024 Koninklijke Philips N.V. All rights are reserved. Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



www.philips.com  
Printed in the Netherlands.

4522 991 88411 \* DEC 2024