

# PHILIPS

## 6800 series

### Sonic electric toothbrush

Smart screen

Advanced sonic technology

Pressure sensor ring

5 brushing modes

HX3792/02



## Multi-dimensional oral care in a smart way

Remove up to 99.9% more plaque\*

Displays everyday brushing performance with build-in smart screen. Gentle on gums with new gum pro mode and advanced pressure sensor ring for oral health. Clean with the G3 Premium Gum Care brush head for an exceptional brushing experience.

### Optimize your brushing

- Smart screen for an instant brushing performance check
- Optimize your brushing with SmarTimer and QuadPacer
- Easy-Start to ease the transition

### Provides superior clean


- Our technology brings you an effective yet gentle clean
- 5 Modes: Clean, Gum Pro, Deep Clean, Sensitive, White
- Removes up to 99.9% more plaque

### Protects your gums

- Pressure sensor ring helps protect your teeth
- G3 brush head to absorb excessive pressure on your gums

## Highlights

### Pressure sensor ring



Your Philips Sonicare is equipped with an advanced sensor that measures the pressure you apply while brushing. If you apply too much pressure, the toothbrush will provide immediate feedback to indicate that you need to adjust the pressure. This feedback is given by a light ring with orange color at the bottom of the handle.

### Advanced Sonic Technology



Our unique sonic sweeping motion is carefully calibrated to create gentle microbubbles that reach deep between the teeth. 31,000 brush strokes per minute\*\* gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

### Smart screen



The built-in smart screen provides instant feedback on your brushing pressure, duration, battery life, and even the notification of when to

replace the brush head. One click to master key messages for better brushing habits without installing an app.

### G3 Premium Gum Care brush head




G3 Premium Gum Care brush head is wrapped with gel. It gently absorbs excessive pressure on your gums and teeth for your best cleaning experience. The top of the bristles is rounded, and the arc design is adopted to increase the contact surface of the bristles, protect the gums and clean them efficiently.

### SmartTimer and QuadPacer



2-minute SmartTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.


### 5 brushing modes



Philips Sonicare 6800 Series has 5 modes for all your brushing needs. Clean mode provides an exceptional daily clean. Gum Pro mode with


triple frequency conversion provides clean and care in one choice. Deep mode provides for invigorating and thorough cleaning. Sensitive mode gently takes care of delicate gums and clean teeth. White mode removes stains efficiently.

### Removes up to 99.9% more plaque



The Sonicare electric toothbrush with advanced sonic technology is clinically proven to remove up to 99.9% more plaque\* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

### Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

# Sonic electric toothbrush

6800 series

HX3792/02

## Specifications

### Power

Voltage: 5Vdc

### Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 14 days

Battery type: Lithium ION

### Design and finishing

Color: Light Pink Metallic

### Service

Warranty: 2-year limited warranty

### Ease of use

Brush head system: Easy click-on brush heads

### Items included

Charger: 1 charger

Handle: 1 6800 Series

Brush heads: 2 G3 Premium Gum Care

Travel Pouch: 1

### Cleaning performance

Timer: BrushPacer and SmarTimer

Pressure feedback: Ring lights up in orange at the bottom of the handle

### Modes

Clean: For exceptional everyday clean (2 minutes)

White: Removes stains on teeth (2 minutes 30 seconds)

Deep Clean: For invigorating and thorough cleaning (3 minutes)

Sensitive: Gentle cleaning for gums (2 minutes)

Gum Pro: Exceptional clean and gentle massage for gums (3 minutes)

### Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

