

PHILIPS

Snoring Relief Band

SmartSleep

86% reported reduced snoring¹

Clinically proven

For positional back snoring

30-day money back guarantee

SN3710/60



Quieter nights for you and your partner¹

The Snoring Relief Band is the first Philips solution designed to alert you before your snoring starts. The sensor is worn in a strap around your torso and delivers gentle vibrations to turn you from your back to your side.

Personalized experience

- Personalized improvements

Easy to read data

- Sensor data on the screen

Easy to experience

- How does it work?
- Vibrating technology
- Fall asleep in a favorite position
- This solution is intended for:

30 day moneyback

- Money back guarantee

Specifications

Product Details

Sensor Dimensions: 3.1 x 1.7 x 0.5 inches

Sensor Weight: 1 oz

Strap Length: 28-51 inches

Strap Weight: 1.16 oz

Material: 42% polyester, 40 % nylon, 18 % spandex

Package contents

Include: Sensor, Chest strap, Charging cable, Instructions for Use

Cleaning

Chest strap: Machine washable, lay flat to air dry

Technical Specifications

Charging: Direct charge, 8 hours

Operation: 8 sleep sessions of use between full charges

Battery: Lithium Polymer 3.7v

