

# PHILIPS

## AVENT



Double Electric  
Breast Pump

Hands-free



SCF547/11

## Free your hands. Free your time.

### Mimics your baby's drinking rhythm

The Philips Avent Hands-free Double Electric Breast Pump is made for effective\* and comfortable expression. It mimics your baby's rhythm and with power, gentleness and lightness, it's your ideal hands-free pump companion.

#### Physical freedom while expressing

- Move freely with the ultra-light collection cups
- Wearable comfort for maximum movement

#### Mimics baby's drinking rhythm

- Express effectively\* with a rhythm that mimics your baby's
- Quiet hospital-strength pumping

#### Easy to use, simple to clean and assemble

- Track milk flow with the fully transparent collection cups
- Easy to clean and assemble with few parts
- Personalise your pumping session

#### Inclusive, seamless fit

- Silicone shields gently shape to your breasts
- A range of breast shield and insert sizes for every mum

## Hands-free

Double Electric Breast Pump

SCF547/11

# Highlights

### Mimics baby's drinking rhythm



Babies know how to drink best! That's why our hands-free electric breast pump mimics their natural drinking rhythm – so you can find your own rhythm for effective\* expression and optimal milk output. Pumping up to 85 times per minute, the Philips Avent Hands-free Breast Pump is 2 x faster than most other pumps\*\*.

### Hospital-strength motor




The nearly silent motor unit delivers the hospital-strength power you expect from a traditional breast pump, yet from a compact and wearable unit. With its rechargeable battery, the cordless motor unit pumps powerfully and discreetly anywhere.

### Unique SkinSense shields



Our unique SkinSense silicone breast shields are designed for ultimate comfort, softly shaping around your breasts using your body's natural warmth. The breast shields are made with food-grade silicone for complete peace of mind.

### Transparent collection



The see-through collection cups provide an unobstructed view inside. They help you correctly place your nipples for your best milk output and visually track the volume of breast milk. For confident pumping from the first drop until finish.

### Ultra-light collection cup



One collection cup weighs just a little over 100 grams when empty - up to 3 x lighter than a fully wearable breast pump. The cups sit comfortably in-bra with barely a sign that they're there.

### Wear it your way




Our pump comes with a braided carrying cord and apparel clip for wearing the motor unit cross-body or attached to your waistband. With the freedom to wear it your way, you stay in your flow knowing control is right by your side. Use the tube coupler and cord adjuster to precisely fine-tune your fit.

### Personalised pumping




Get the most from every session with 8 suction levels for stimulating the flow of milk and 16 suction levels for expressing it. You can fine-tune as you pump and with your ideal setting automatically memorised, your pump is always ready to go.

### Easy to clean and assemble



We designed our pump and accessories to be as simple as possible. With only a few connecting pieces, each collection cup is a breeze to take apart, clean (by hand or in a dishwasher), disinfect and put back together.

### Find the right fit



We support every mum to feel their best and having the right fit is an essential part of that. Our range of breast shields and shield insert sizes fits 99% of mums. You can experience your ultimate comfort while pumping. Available to buy separately.

# Hands-free

Double Electric Breast Pump

SCF547/11

## Specifications

### Power

**Source:** Rechargeable battery / USB-C (adapter not included)

### Material

**Cup, front cap, breast shield cover:**

Polypropylene

**Breast shield, membrane, valve, insert:** Liquid silicone rubber

**Silicone tube:** Silicone

### What is included

**Breast milk bags:** 25 pcs

**Carrying cord:** 1 pcs

**Tube coupler:** 1 pcs

**Rechargeable motor unit (with a detachable**

**clip):** 1 pcs

**Collection cup (assembled):** 2 pcs

**Silicone tubes:** 2 pcs

**Breast shield (dust) cover:** 2 pcs

**Breast shield (21 mm, 24 mm):** 2 pcs each

**125 ml (4 oz) Natural Response bottle (with**

**Flow 2 teat):** 1 pcs

**Nipple sizing card:** 1 pcs

**Insert 19 mm:** 2 pcs

**USB cable:** 1 pcs

### Ease of use

**Breast pump use:** Easy to use, simple to clean and assemble

### Functions

**Flexibility in movement:** Hands-free pumping

**Personalisation settings:** Stimulation mode with

8 suction levels and expression mode with

16 suction levels

**Motor unit:** Rechargeable battery / USB-C, Fully charging the battery takes up to 1.5 hours and lasts approx. 5 pumping sessions

**Breast shield and insert specifics:** Several sizes available for the best fit

